

## General Description

The MAX32674C may be programmed as an algorithm hub (AlgoHub) or sensor hub (SensorHub) which provides the following innovative features:

- ▶ Blood Pressure Trending (SensorHub) Heart rate, and SpO2 algorithms developed by a team of machine learning, data scientist, algorithm experts.
- ▶ AFE sensor drivers.
- ▶ OTA update compatibility which allows the flashing of the latest algorithm .msbl file.
- ▶ Low-powered SensorHub to manage AFE and accelerometer data collection.
- ▶ Faster time to market; development time cut by at least six months.
- ▶ Sample host code reduces integration time.
- ▶ Reference design includes a host processor, enclosure and Android GUI app which allows for wrist-band evaluation.

## Evaluation Platform

The MAXREFDE104# may be used to evaluate biometric algorithms. The MAXREFDES104# platform includes the following major components

- ▶ MRD104 Micro Board
  - 32MB Flash Memory (MX25U51245GZI54).
  - MAX32670 Microcontroller pre-programmed with MAX32674C bootloader firmware. (SensorHub when flashed with .msbl file)
  - MAX20360 Power Management IC (PMIC)
  - MAX32666 Host Microcontroller (Host MCU) with Integrated BLE
  - Status LEDs

- ▶ MAX86176 Sensor Board
  - MAX86176 ECG and PPG Biosensor AFE
  - LEDs: 1 Green, 1 Red, 1 IR [OSRAM SFH7016]
  - 3 Photodiodes [Vishay VEMD8080]
  - 3-Axis Accelerometer [ST Micro LIS2DS12]
- ▶ Li-Po Battery 105mAh [Adafruit 1570]
- ▶ MAXDAP-TYPE-C Programming Adapter
- ▶ Cy Smart CY5677 BLE USB Dongle
- ▶ Stainless Steel (ANSI SS304) ECG Dry Electrodes (3 pcs.).
- ▶ 3D Printed Enclosure and Band



Figure 1. MAXREFDES104#

## Contents

General Description .....	1	I2C SensorHub Output FIFO Format ...	58
Evaluation Platform.....	1	SensorHub, Heart Rate, SpO2, AEC/AGC Report Configuration.....	62
MAXREFDES104# Architecture.....	3	SensorHub, PPG Raw Report Configuration .....	65
AlgoHub, SensorHub Configurations .....	3	SensorHub BPT, Heart Rate, SpO2 Report Configuration.....	66
AlgoHub Configuration .....	3	BPT User Calibration Vector Report .	67
SensorHub Configuration .....	4	BPT Estimation Report and Heart Rate, SpO2 Report .....	69
When to use AlgoHub configuration or SensorHub configuration.....	5	BPT Restore User Calibration Vectors .....	72
MAX32674C Reset .....	6	Configuring a New MAX32674C.....	72
Reset to Bootloader or Application Mode	6	Application .msbl Programming Sequence 73	
Reset to an Unsupported State .....	7	SpO2 Coefficients for Final Factor.....	82
AlgoHubSensorHub Handshaking: I2C and MFIO pin .....	8	APIs for Sleep, Shutdown.....	82
AlgoHubSensorHub I2C Communications..	8	MAX32674C Processing Capabilities .....	83
MAX32674C I2C Bit Transfer Process....	9	Heart Rate Algorithm Performance.....	83
Read Status Byte .....	10	SpO2 Algorithm Performance .....	84
SensorHub I2C Command Response Definitions .....	10	References.....	85
Configuring the AlgoHub and Reading Processed Data .....	44		
AlgoHub Output FIFO Format.....	44		
AlgoHub AEC Normal Algorithm Report Configuration and Batched Processing.	46		
AlgoHub AGC Normal Algorithm Report Configuration .....	52		
Configuration for AlgoHub Normal Algorithm Report when Host Implements AEC/AGC .....	55		
AlgoHub Processing Response Duration .....	58		
Configuring the SensorHub and Reading Processed Data .....	58		

## MAXREFDES104# Architecture

The MAX32674C (AlgoHubSensorHub) is pre-programmed with a bootloader that accepts Analog Devices biometric algorithms, and the bootloader is used to flash the algorithm to its application flash memory. The algorithms are packaged as .msbl files. The MAX32674C is not factory pre-programmed with the application algorithm .msbl. Please see the MAXREFDES104# User Guide for instructions on how to flash the .msbl MAX32674C (The initial stock of MAXREFDES104 comes populated with MAX32670 which has been programmed with the bootloader so that it behaves as a MAX32674C).

The MAXREFDES104# includes the MAX32666 which acts as the host processor. The host processor is responsible for initializing AlgoHubSensorHub and reading the biometric data from the AlgoHubSensorHub FIFO. The host processor on the MAXREFDES104# is also programmed to send BLE data to an Android app. Follow the instructions MAXREFDES104# User Guide to program the .bin to the MAX32666.

**When using the MAXREFDES104, the following three software and firmware must all be updated using the latest and same software release .zip package. Please refer to the MAXREFDES104 User Guide for detailed software, firmware installation procedure.**

1. Flash the micro board **.bin** file using drag & drop to the DAPLINK folder.
2. Flash the algorithm **.msbl** file to MAX32674C using the batch file.
3. Install the MAXREFDES104 PC GUI **.exe** file. (if evaluating BPT, install Wellness app **.apk** file to an Android device and grant all permissions.)

### AlgoHub, SensorHub Configurations

The hardware and firmware for the MAX32674C may be configured in the AlgoHub or SensorHub configurations. Both configurations require accel data

#### AlgoHub Configuration

In the AlgoHub configuration, the MAX32666 host is connected to the MAX86176 AFE and the LIS2DS accelerometer. The host drivers are responsible for the setting of the AFE, accel registers and for retrieving the PPG, accel data. After the MAX32674C has been initialized to AlgoHub mode via I2C commands, the host pipes in the PPG, accel data to the MAX32674C algorithms and the MAX32674C provides heart rate and SpO2 data to the host. The MAX32674C is power optimized to deep sleep when it is idle. By using the combined ECG, PPG MAX86176 IC, the host processor can retrieve time synchronized ECG and PPG data.




Figure 2. AlgoHub Architecture diagram for MAXREFDES104# health-sensing platform

### SensorHub Configuration

In the SensorHub configuration, the MAX32674C is connected to the MAX86176 AFE and the LIS2DS12 accelerometer. After the MAX32674C has been initialized to SensorHub mode via I2C commands, then the MAX32674C drivers will handle the AFE, accel register settings and retrieve the PPG, accel data. The MAX32674C algorithms process the raw data to provide BPT, heart rate and SpO2 data to the host. The MAX32674C is power optimized to deep sleep when it is idle.




Figure 3. SensorHub Architecture diagram for MAXREFDES104# health-sensing platform.

### When to use AlgoHub configuration or SensorHub configuration

The AlgoHub configuration should be used when the system requirements require that ECG and PPG run simultaneously or when it is required to have accel, ECG and PPG data time synchronized. The AlgoHub mode requires more host integration effort for AFE, accel register settings and for piping in the PPG and accel data to the MAX32674C and for creating the host code which responds to the AlgoHub algorithm requests to update the PPG register settings of the AFE. The total system power will be more for AlgoHub configuration since there will be more host to SensorHub communications and the host must handle the dynamic management of the PPG AFE and accel sensors.

The SensorHub configuration should be used when the system requirements do not require that ECG and PPG run simultaneously – ECG and PPG are run in stand-alone modes. The SensorHub mode requires less host integration effort since the MAX32674C includes the drivers to retrieve the PPG, accel data from the MAX86175 and LIS2DS12. The total system power will be less for SensorHub configuration since there will be less host to SensorHub communications and the dynamic management of the PPG AFE and the accel is off-loaded to the power optimized MAX32674C. ECG is currently not demonstrated in SensorHub configuration.


## MAX32674C Reset

### Reset to Bootloader or Application Mode

Two GPIO pins from the host are needed to control the RSTN and multifunction input/output (MFIO) pins.

To enter Bootloader mode:

- ▶ Set the RSTN pin low.
- ▶ While RSTN is low, set the MFIO pin to low and set either SWDCLK or UART0\_RX to high. (Do this at least 1ms before the RSTN pin is set to high.)
- ▶ After the 10ms has elapsed, set the RSTN pin to high.
- ▶ After an additional 50ms has elapsed, the SensorHub is in Bootloader mode.



*Figure 4 Entering bootloader mode using the RSTN pin and the MFIO GPIO pin.*

To enter Application mode (The MAX32674C is not pre-programmed with the algorithm application .msbl):

- ▶ Set the RSTN pin low.
- ▶ While RSTN is low, set the MFIO pin to high and set either SWDCLK or UART0\_RX to high. (Do this at least 1ms before the RSTN pin is set to high.).
- ▶ After the 10ms has elapsed, set the RSTN pin to high



- ▶ After an additional 50ms has elapsed, the SensorHub is in Application mode and the application performs its initialization of the application software.
- ▶ After approximately 1.5 second from when the RSTN pin was set to high, the application completes the initialization, and the device is ready to accept I<sup>2</sup>C commands.
- ▶ Set the ACCEL\_CS\_N pin to low.




Figure 5 Entering application mode using the RSTN pin and MFIO pin.

To enter Application mode by timing out from Bootloader mode:

- ▶ Set the RSTN pin low.
- ▶ While RSTN is low, set the MFIO pin to low and set either SWDCLK or UART0\_RX to high. (Do this at least 1ms before the RSTN pin is set to high.)
- ▶ After the 10ms has elapsed, set the RSTN pin to high.
- ▶ After an additional 50ms has elapsed, the AlgoHubSensorHub is in Bootloader mode.
- ▶ If no I<sup>2</sup>C commands are sent to the AlgoHubSensorHub within the next 1s, then the AlgoHubSensorHub will automatically switch to application mode.




Figure 6 Entering application mode by timing out from Bootloader mode.

## Reset to an Unsupported State

Do not set P0.1 (ACC\_INT) and P0.8 (ACC\_CS) both low when powering up or exiting reset – doing so will cause the MAX32674C to enter an unsupported state. Set either or both P0.1, P0.8 to high when powering up or exiting reset.

## AlgoHubSensorHub Handshaking: I2C and MFIO pin

Normally, when the AlgoHubSensorHub is idle, it switches to deep sleep mode to save power. An external interrupt-like sensor, host MFIO, or RTC alarm forces the AlgoHubSensorHub to wake up.

The host is required to wake up the AlgoHubSensorHub prior to any I2C communication by:

- ▶ Setting the MFIO pin to low at least 300µs before the beginning of an I2C transaction to wake the AlgoHubSensorHub.
- ▶ Keeping the MFIO pin low during the I2C transaction.
- ▶ Setting MFIO to high after the end of I2C communication to allow the AlgoHubSensorHub to switch back to deep sleep.




Figure 7. I2C and MFIO pin: host uses MFIO for enabling host communications.

## AlgoHubSensorHub I2C Communications

A host uses the I2C bus to communicate with the AlgoHubSensorHub (slave) using a series of commands. The default CMD\_DELAY is 2ms. The I2C data rate 400Kbps has been tested; the 1000Kbps data rate has not been tested.

Slave\_WriteAddress and Slave\_ReadAddress are set to 0xAA and 0xAB, respectively.

A generic write command includes the following fields:

```
Slave_WriteAddress(1 byte) | Command_Family(1 byte) | [Index byte] | [Write
byte] | [additional command byte(s)]
```

A generic response includes the following fields:

```
Slave_ReadAddress(1 byte) | Read Status Byte | Value (multiple bytes)
```






Figure 8. Example I2C, MFIO example command response.

### MAX32674C I2C Bit Transfer Process.

The read status byte is an indicator of success (0x00) or failure, as detailed in the following table. The defined bit transfer process is described below. It is recommended that I2C GPIO 'bit-bang' software be implemented on the host if the host MCU I2C hardware/HAL is not compatible with AlgoHubSensorHub protocol.

Both SDA and SCL signals are open-drain circuits. Each has an external pullup resistor that ensures each circuit is high when idle. The I2C specification states that during data transfer, the SDA line can change state only when SCL is low, and that SDA is stable and able to be read when SCL is high. Typical I2C write/read transactions are shown below.




Figure 1. I2C Write/Read data transfer from host microcontroller.

The read status byte is an indicator of the success or failure of the Write Transaction. The read status byte must be accessed after each write transaction to the device. This ensures that write transaction processing is understood and any errors in the device command handling can be corrected. The value of the read status byte is summarized in the next table.

### Read Status Byte

The read status byte is an indicator of success (0x00) or failure, as detailed in the following table.

**Table 1. Read Status Byte Definition**

READ STATUS BYTE VALUE	DESCRIPTION
0x00	Application mode: The write transaction was successful.
0x01	Illegal Family Byte and/or Index Byte was used. Verify that the Family Byte, Index Byte are valid for the host command sent.
0x02	Illegal Index Byte and/or Write Byte was used. Verify that the Index Byte and Write Byte(s) are valid for the host command sent.
0x03	Incorrect number of bytes sent for the requested Family Byte. Verify that the correct number of bytes are sent for the host command.
0x04	Illegal configuration value was attempted to be set. Verify that the Family Byte, Index Byte, and Write Byte are correct.
0x05	Not used in application mode. (In bootloader: Device is busy. Insert delay and resend the host command.)
0x80	Not used. General error while receiving/flashing a page during the bootloader sequence. Not used.
0x81	Bootloader checksum error while decrypting/checking page data. Verify that the keyed .msbl file is compatible with MAX32674C.
0x82	Bootloader authorization error. Verify that the keyed .msbl file is compatible with MAX32674C.
0x83	Bootloader detected that the application is not valid.
0x84	Bootloader: application program was not erased before trying to flash a new one.
0xAA	Bootloader mode: The write transaction was successful.
0xAB	Bootloader: partial page data received; there is still more page data.
0xFE	Device is busy. Try again. Increase the delay before the command and increase the CMD_DELAY.
0xFF	Unknown error. The AlgoHubSensorHub is in deep sleep unless the host sets the MFIO pin low 300us before during the I2C communications. When switching to bootloader mode, allow 50ms for initialization. When switching to application mode, allow 1.5s for initialization.
NAK	NAK received. AlgoHubSensorHub was busy. Resend command after 1ms with a maximum of five retries. If this issue persists, then empty the FIFO by reading all the data. Verify that the hardware I2C/MFIO rise times, voltage levels, and grounding are correct. Verify that the MFIO line fall time is clean; increase the MFIO pin low time to wake to 300us

### SensorHub I2C Command Response Definitions

The following table defines the host to MAX32647C I2C command response definitions. The Read Status Byte is not listed in this table.

**Table 2. AlgoHubSensorHub I2C Command Response Definitions**

FAMILY BYTE	INDEX BYTE	WRITE BYTE(S)	DESCRIPTION	READ BYTE(S) RESPONSE
0x00	0x00	-	Read sensor AlgoHubSensorHub status.	<b>Err0[0]</b> : 0 = No error; 1 = Sensor communication problem <b>Err1[0]</b> : Not used <b>Err2[0]</b> : Not used <b>DataRdyInt[3]</b> : 0 = FIFO below threshold; 1 = FIFO filled to threshold or above. <b>FifoOutOvrInt[4]</b> : 0 = No FIFO overflow; 1 = AlgoHubSensorHub output FIFO overflowed, data lost. <b>FifoInOvrInt[5]</b> : 0 = No FIFO overflow; 1 = AlgoHubSensorHub Input FIFO overflowed, data lost.
0x01	0x00	<b>0x00</b> : Exit bootloader mode, enter application mode. (CMD_DELAY = 1.5s) <b>0x01</b> : Shutdown the MAX32674C. Restart by power cycling or pulsing RSTN (application mode) <b>0x08</b> : Enter bootloader mode (if commanded within the first 20ms after reset)	Set the device operating mode.	-
0x02	0x00	-	Read the device operating mode.	0x00: Application operating mode. 0x02: Reset.

				0x08: Bootloader operating mode.
0x10	0x00	<b>0x00:</b> (no data) <b>0x01:</b> Sensor Data <b>0x02:</b> Algorithm Data <b>0x03:</b> Sensor Data and Algorithm Data <b>0x04:</b> (no data) <b>0x05:</b> Sample Counter byte, Sensor Data <b>0x06:</b> Sample Counter byte, Algorithm Data <b>0x07:</b> Sample Counter byte, Sensor Data and Algorithm Data	Set the output format of the AlgoHubSensorHub.	
0x10	0x01	<b>0x01 to 0xFF:</b> AlgoHubSensorHub Interrupt Threshold for FIFO.	Set the threshold for the FIFO interrupt bit. The bit DataRdyInt of the AlgoHubSensorHub status byte is set when this threshold is reached.	-
0x11	0x01	-	Read the threshold for the FIFO interrupt bit. The bit DataRdyInt of the AlgoHubSensorHub status byte is set when this threshold is reached.	<b>0x01 to 0xFF:</b> AlgoHubSensorHub Interrupt Threshold for FIFO.
0x10	0x02	<b>0x01 to 0xFF:</b> Isb is 40ms. N, where a samples report is generated once every N samples.  BPT Isb is 8ms.	Set the samples report period (e.g., a value of 25 means a samples report is generated once every 25 samples).  BPT, the minimum samples report period is 4. (e.g., a value of 5 means that the report period is 40ms.	-
0x11	0x02	-	Read the samples report period (e.g., a value of 25 means a samples report is generated once every 25 samples).	<b>0x01</b> (default) to <b>0xFF:</b> Isb is 40ms. N, where a samples report is generated once every N samples.

0x10	0x03	<b>0x02 to 0xFF:</b> New I <sup>2</sup> C address (8-bit I <sup>2</sup> C write address)	Change I <sup>2</sup> C address of the MAX32674C	-
0x10	0x04	<b>0x00 to 0xFF:</b> Counter	Set the AlgoHubSensorHub counter.	-
0x11	0x04	-	Read the AlgoHubSensorHub counter.	0x00 to <b>0xFF</b> : Counter
0x11	0x05	-	Read PPG output FIFO samples report size	Number of bytes in the PPG samples report
0x12	0x00	-	Get the number of samples available in the output FIFO	Number of samples available in the FIFO.
0x12	0x01 (CMD_DELAY = 5ms)	-	Read data stored in output FIFO.	Samples Report from Output FIFO. The internal FIFO read pointer increments once the sample size bytes have been read. See Samples Report Table for more details.
0x14 AlgoHub	0x00	AlgoHub Input Frame 1 AlgoHub Input Frame 2 . . . AlgoHub Input Frame N-1 AlgoHub Input Frame N  <b>N : [1 to 25]</b>  CMD_DELAY = 16ms	Write AlgoHub input frame data to the AlgoHub input FIFO.  <b>AlgoHub Input Frame (24 Bytes) Content:</b>  PPG1: 3 byte PPG value: [MSB .. LSB] for WHRM HR Channel 1  PPG2: 3 byte PPG value: [MSB .. LSB] for WHRM HR Channel 2  PPG3: 3 byte PPG value: [MSB .. LSB]for SpO2 IR Channel  PPG4 : 3 byte PPG value: [MSB .. LSB] for SpO2 Red Channel  PPG5 : 3 byte PPG value:[MSB .. LSB] for N/A  PPG6 : 3 byte PPG value: [MSB .. LSB] for N/A  Accelerometer set to ±8g:  ACCLX: 2-byte accel X value (2's complement .001g): [MSB LSB]  ACCLY: 2-byte accel Y value (2's complement .001g): [MSB LSB]	mm nn  mmnn is the ppg, accel number of bytes received (16-bit unsigned)s

			ACCLZ: 2-byte accel Z value (2's complement .001g): [MSB LSB]	
0x42 Sensor Hub	0x04		Read the attributes of the accelerometer.	Number of bytes in a word for this sensor, Number of registers available for this sensor.
0x42 Sensor Hub	0x06		Read the attributes of the MAX86176.	Number of bytes in a word for this sensor, Number of registers available for this sensor.
0x40 Sensor Hub	0x04	[reg_addr] [reg value]	Write a value to a writable accelerometer sensor register.	
0x41 Sensor Hub	0x04	[reg_addr]	Read accelerometer sensor register.	Register value (byte)
0x40 Sensor Hub	0x06	[reg_addr] [reg value]	Write a value to a writable MAX86176 register.	
0x41 Sensor Hub	0x06	[reg_addr]	Read the value of a MAX86176 register.	Register value (byte)
0x44 Sensor Hub	0x04	<b>0x00, 0x00:</b> Disable sensor hub accelerometer  CMD_DELAY = 50ms <b>0x01, 0x00:</b> Enable sensor hub accelerometer CMD_DELAY = 50ms	Write the enable/disable the accelerometer.	-
0x44 Sensor Hub	0x06	<b>0x00:</b> Disable, CMD_DELAY = 200ms <b>0x01, 0x00:</b> Enable Raw Data, CMD_DELAY = 500ms	Write the enable/disable the MAX86176 AFE.	-
0x44 AlgoHub	0x07	<b>0x00, 0x01:</b> Disable WAS AlgoHub, external sensor input, CMD_DELAY = 200ms. <b>0x01, 0x01:</b> Enable WAS	Enable/Disable Wearable Algorithm Suite AlgoHub.	-



		AlgoHub, enable external sensor input . CMD_DELAY = 500ms		
0x46 for write. AlgoHub	0x07	0x00 [A_MSB ... A_LSB] [B_MSB ... B_LSB] [C_MSB ... C_LSB]	Write SpO2 calibration coefficients (12 bytes comprised of three 32-bit signed values, scaled up by 100,000)	-
0x47 for read. AlgoHub	0x07	0x00	Read SpO2 calibration coefficients (12 bytes comprised of three 32-bit signed values, scaled up by 100,000) Default: 0xFFFF9482C FFE3792B 00A8E9A4 A = -4.4027606 (0xFFFF9482C) B = -18.695252 (0xFFE3792 ) C = 110.6986082 (0x00A8E9A4)	32-bit signed integer A, 32-bit signed integer B, 32-bit signed integer C  Values scaled up by 100,000
0x46 for write. AlgoHub	0x07	0x01 [MSB LSB]	Write SpO2 motion-detection period (unsigned 16-bit int, seconds). The WAS algorithm will consider the state to be motionless if the motion is below the threshold for this duration of time.	-
0x47 for read. AlgoHub	0x07	0x01	Read SpO2 motion-detection period. The WAS algorithm will consider the state to be motionless if the motion is below the threshold for this duration of time. Default: 0x0002	MSB of period, LSB of period (16-bit unsigned integer, seconds)
0x46 for write. AlgoHub	0x07	0x02 [MSB ... LSB]	Write SpO2 motion-detection threshold (signed 32-bit int, equal to 10 <sup>5</sup> x milli-g threshold value)	-
0x47 for read. AlgoHub	0x07	0x02	Read SpO2 motion-detection threshold (signed 32-bit int, equal to 10 <sup>5</sup> x milli-g threshold value) Default: 0x01C9C380 (0.3g)	4 bytes (32-bit signed integers which are the milli-g motion threshold times 100,000)
0x46 for write. AlgoHub	0x07	0x03 [byte]	Write SpO2 AGC Timeout (sec) .	
0x47 for read. AlgoHub	0x07	0x03	Read SpO2 AGC Timeout (sec) . Default: 0x1E	SpO2 AGC Timeout (8-bit unsigned, seconds)
0x46 for write. AlgoHub	0x07	0x04	Write the timeout duration for SpO2 measurement in seconds (1 byte).	-

0x47 for read. AlgoHub	0x07	0x04	Read the timeout duration for SpO2 measurement in seconds. Default: 0x78	SpO2 algorithm timeout (8-bit unsigned, seconds)
0x46 for write. AlgoHub	0x07	0x05 [byte]	Write initial HR algorithm value (8-bit unsigned).	-
0x47 for read. AlgoHub	0x07	0x05	Read initial HR algorithm value. Default: 0x3C	Initial heart rate (8-bit unsigned)
0x46 for write. AlgoHub	0x07	0x06 [MSB] [LSB]	Write height (16-bit unsigned, cm).	-
0x47 for read. AlgoHub	0x07	0x06	Read height. Default: 0x00AF	Height (16-bit unsigned, cm)
0x46 for write. AlgoHub	0x07	0x07 [MSB] [LSB]	Write weight (16-bit unsigned, kg).	-
0x47 for read. AlgoHub	0x07	0x07	Read weight. Default: 0x004E	Weight (16-bit unsigned, kg)
0x46 for write. AlgoHub	0x07	0x08 [byte]	Write age (8-bit unsigned, years).	
0x47 for read. AlgoHub	0x07	0x08	Read age. Default: 0x1E	Age (8-bit unsigned, years)
0x46 for write. AlgoHub	0x07	0x09 [byte]	Write gender. Byte: <b>0x00</b> : Male <b>0x01</b> : Female	
0x47 for read. AlgoHub	0x07	0x09	Read gender	Gender <b>0x00</b> : Male (default) <b>0x01</b> : Female
0x46 for write. AlgoHub	0x07	0x0A [WAS algo operation mode byte]	Set the WAS algorithm operation mode (can be switched in runtime): <b>0x00</b> : Continuous HRM + Continuous SpO2 (default). <b>0x01</b> : Continuous HRM + One-Shot SpO2	-

			<b>0x02:</b> Continuous HRM <b>0x03:</b> Sampled HRM <b>0x04:</b> Sampled HRM + One-Shot SpO <sub>2</sub> <b>0x05:</b> Activity Tracking ONLY <b>0x06:</b> SpO <sub>2</sub> Calibration Data Collection	
0x47 for read. AlgoHub	0x07	0x0A	Read the WAS algorithm operation mode.	<b>0x00:</b> Continuous HRM, continuous SpO <sub>2</sub> (default) <b>0x01:</b> Continuous HRM, one-shot SpO <sub>2</sub> <b>0x02:</b> Continuous HRM <b>0x03:</b> Sampled HRM <b>0x04:</b> Sampled HRM, one-shot SpO <sub>2</sub> <b>0x05:</b> Activity tracking only <b>0x06:</b> SpO <sub>2</sub> calibration
0x46 for write. AlgoHub	0x07	0x0B [byte]	Write the enable/disable AEC byte: <b>0x00:</b> Disable <b>0x01:</b> Enable	-
0x47 for read. AlgoHub	0x07	0x0B	Read the enable/disable AEC.	<b>0x00:</b> Disabled <b>0x01:</b> Enabled (default)
0x46 for write. AlgoHub	0x07	0x0C [byte]	Write the Skin Contact Detection (SCD) algorithm enable <b>0x00:</b> Disable <b>0x01:</b> Enable	-
0x47 for read. AlgoHub	0x07	0x0C	Read the Skin Contact Detection (SCD) algorithm enable	<b>0x00:</b> Disable <b>0x01:</b> Enable (default)
0x46 for write. AlgoHub	0x07	0x0D [MSB] [LSB]	Write adjusted target PD current period (16-bit unsigned, seconds )	-
0x47 for read. AlgoHub	0x07	0x0D	Read adjusted target PD current period in seconds. Default: 0x0708	Adjusted target PD current period (16-bit unsigned, seconds)
0x46 for write. AlgoHub	0x07	0x0E [byte]	Write HR motion magnitude threshold (16-bit unsigned, 0.001g)	-
0x47 for read. AlgoHub	0x07	0x0E [MSB] [LSB]	Read HR motion magnitude threshold. Default: 0x0032 (0.05g)	[MSB] [LSB] motion magnitude threshold (16-bit unsigned, 0.001g)
0x46 for write. AlgoHub	0x07	0x0F [MSB] [LSB]	Write minimum PD current (16-bit unsigned, 0.1uA).	-

0x47 for read. AlgoHub	0x07	0x0F	Read minimum PD current. Default: 0x0032	Minimum PD current (16-bit unsigned, 0.1uA)
0x46 for write. AlgoHub	0x07	0x10 [MSB] [LSB]	Write initial PD current (16-bit unsigned, 0.1uA). This sets the target PD current you would like AEC algorithm to maintain initially. It does not correspond to any register. Once you set what PD current you need, algorithm will calculate the appropriate LED current.	-
0x47 for read. AlgoHub	0x07	0x10	Read initial PD current. Default: 0x0064	Initial PD current (16-bit unsigned, 0.1uA)
0x46 for write. AlgoHub	0x07	0x11 [MSB] [LSB]	Write target PD current (16-bit unsigned, 0.1uA). Applicable only if Auto Target PD Current Calculation is enabled.	
0x47 for read. AlgoHub	0x07	0x11	Read target PD current. Default: 0x0064	Target PD current (16-bit unsigned, 0.1uA)
0x46 for write. AlgoHub	0x07	0x12 [byte]	Write enable/disable automatic calculation of target PD current. byte: <b>0x00</b> : Disable <b>0x01</b> : Enable	-
0x47 for read. AlgoHub	0x07	0x12	Read enable/disable automatic calculation of target PD current	<b>0x00</b> : Disable <b>0x01</b> : Enable (default)
0x46 for write. AlgoHub	0x07	0x13 [byte1 byte 2]	Write minimum AFE integration time. byte1: <b>0x00</b> to <b>0x08</b> : PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A) byte2: <b>0x00</b> : 14.8μs <b>0x01</b> : 29.4μs <b>0x02</b> : 58.7μs <b>0x03</b> : 117.3μs	-
0x47 for read. AlgoHub	0x07	0x13	Read minimum AFE integration time.	Minimum AFE integration time for each PPG measurement 1 to 9. [byte1 to byte9] byteN: <b>0x00</b> : 14.8μs (default) <b>0x01</b> : 29.4μs <b>0x02</b> : 58.7μs

				<b>0x03: 117.3μs</b>
0x46 for write. AlgoHub	0x07	0x14 [byte1 byte2]	Write minimum sampling. byte1: <b>0x00 to 0x08</b> : PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A)  byte2: <b>0x01 to 0x80</b>	-
0x47 for read. AlgoHub	0x07	0x14	Read minimum sampling rate and averaging.	Minimum sampling rate for each PPG measurement 1 to 9. [byte1 to byte9]  byteN: <b>0x01 to 0x80</b> (default 0x01)
0x46 for write. AlgoHub	0x07	0x15 [byte1 byte2]	Write maximum integration time: <b>0x00 to 0x08</b> : PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A)  byte2: <b>0x00</b> : 14.8μs <b>0x01</b> : 29.4μs <b>0x02</b> : 58.7μs <b>0x03</b> : 117.3μs	-
0x47 for read. AlgoHub	0x07	0x15	Read maximum integration time:	Maximum integration time for each PPG measurement 1 to 9. [byte1 to byte9]  byteN: <b>0x00</b> : 14.8μs <b>0x01</b> : 29.4μs <b>0x02</b> : 58.7μs <b>0x03</b> : 117.3μs (default)
0x46 for write. AlgoHub	0x07	0x16 [byte1 byte 2]	Write maximum sampling rate and averaging: byte1: <b>0x00 to 0x08</b> : PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A)  byte2: <b>0x01 to 0x80</b>	-
0x47 for read. AlgoHub	0x07	0x16	Read maximum sampling rate.	Maximum sampling rate for each PPG measurement 1 to 9. [byte1 to byte9]  byteN: <b>0x01 to 0x80</b> (default 0x04)
0x46 for write. AlgoHub	0x07	0x17 0xWX, 0xYZ	Write slot and PD configuration for the two HR inputs to the WHRM algorithm.	-

			<p>0xWX, 0xYZ:</p> <p>WX is input 1 of the WHRM algorithm.  W (MS nibble of WX):  W = 0 for Slot 1  W = 1 for Slot 2  W = 2 for Slot 3  W = 3 for Slot 4  W = 4 for Slot 5  W = 5 for Slot 6  W = 7 for Slot not used  X (LS nibble of WX):  X = 0 for PD1  X = 1 for PD2  X = 3 for PD not used.</p> <p>YZ is input 2 of the WHRM algorithm.  Y (MS nibble of YZ):  Y = 0 for Slot 1  Y = 1 for Slot 2  Y = 2 for Slot 3  Y = 3 for Slot 4  Y = 4 for Slot 5  Y = 5 for Slot 6  Y = 7 for Slot not used  Z (LS nibble of YZ):  Z = 0 for PD1  Z = 1 for PD2  Z = 3 for PD not used.</p>	
0x47 for read. AlgoHub	0x07	0x17	<p>Read Slot and PD configuration for the two HR inputs to the WHRM algorithm.</p>	<p>0xWX, 0xYZ  Default: 0x0001</p> <p>WX is input 1 of the WHRM algorithm.  W (MS nibble of WX):  W = 0 for Slot 1  W = 1 for Slot 2  W = 2 for Slot 3  W = 3 for Slot 4  W = 4 for Slot 5  W = 5 for Slot 6  W = 7 for Slot not used  X (LS nibble of WX):  X = 0 for PD1  X = 1 for PD2  X = 3 for PD not used.</p> <p>YZ is input 2 of the WHRM algorithm.  Y (MS nibble of YZ):  Y = 0 for Slot 1  Y = 1 for Slot 2  Y = 2 for Slot 3  Y = 3 for Slot 4</p>



				Y = 4 for Slot 5 Y = 5 for Slot 6 Y = 7 for Slot not used Z (LS nibble of YZ): Z = 0 for PD1 Z = 1 for PD2 Z = 3 for PD not used.
0x46 for write. AlgoHub	0x07	0x18 0xWX 0xYZ	Write Slot and PD configuration for the IR, red inputs to the WSpO2 algorithm.  0xWX, 0xYZ:  WX is the LED/PD used for IR for the WSpO2 algorithm. W (MS nibble of WX): W = 0 for Slot 1 W = 1 for Slot 2 W = 2 for Slot 3 W = 3 for Slot 4 W = 4 for Slot 5 W = 5 for Slot 6 W = 7 for Slot not used X (LS nibble of WX): X = 0 for PD1 X = 1 for PD2 X = 3 for PD not used.  YZ is the LED/PD used for red for the WSpO2 algorithm. Y (MS nibble of YZ): Y = 0 for Slot 1 Y = 1 for Slot 2 Y = 2 for Slot 3 Y = 3 for Slot 4 Y = 4 for Slot 5 Y = 5 for Slot 6 Y = 7 for Slot not used Z (LS nibble of YZ): Z = 0 for PD1 Z = 1 for PD2 Z = 3 for PD not used	-
0x47 for read. AlgoHub	0x07	0x18	Read Slot and PD configuration for the IR, red inputs to the WSpO2 algorithm.	0xWX, 0xYZ Default: 0x1020  WX is the LED/PD used for IR for the WSpO2 algorithm. W (MS nibble of WX): W = 0 for Slot 1 W = 1 for Slot 2 W = 2 for Slot 3 W = 3 for Slot 4 W = 4 for Slot 5 W = 5 for Slot 6

				<p>W = 7 for Slot not used  X (LS nibble of WX):  X = 0 for PD1  X = 1 for PD2  X = 3 for PD not used.</p> <p>YZ is the LED/PD used for red for the WSpO<sub>2</sub> algorithm.  Y (MS nibble of YZ):  Y = 0 for Slot 1  Y = 1 for Slot 2  Y = 2 for Slot 3  Y = 3 for Slot 4  Y = 4 for Slot 5  Y = 5 for Slot 6  Y = 7 for Slot not used  Z (LS nibble of YZ):  Z = 0 for PD1  Z = 1 for PD2  Z = 3 for PD not used</p>
0x46 for write. AlgoHub	0x07	0x1A [byte1 byte2]	<p>Write initial AFE integration time.  byte1:  <b>0x00</b> to <b>0x08</b>: PPG measurement 1 to 9.  (1: Green, 2: IR, 3: red, N/A, ... N/A)</p> <p>byte2:  <b>0x00</b>: 14.8μs  <b>0x01</b>: 29.4μs  <b>0x02</b>: 58.7μs  <b>0x03</b>: 117.3μs</p>	-
0x47 for read. AlgoHub	0x07	0x1A	<p>Read initial AFE integration time.</p>	<p>Initial AFE integration time for each PPG measurement 1 to 9.  [byte1 to byte9]</p> <p>byteN:  <b>0x00</b>: 14.8μs  <b>0x01</b>: 29.4μs  <b>0x02</b>: 58.7μs  <b>0x03</b>: 117.3μs (default)</p>
0x46 for write. AlgoHub	0x07	0x1B [byte1 byte2]	<p>Write initial sampling rate.  <b>0x00</b> to <b>0x08</b>: PPG measurement 1 to 9.  (1: Green, 2: IR, 3: red, N/A, ... N/A)</p> <p>byte2:  <b>0x01</b> to <b>0x80</b></p>	
0x47 for read. AlgoHub	0x07	0x1B	<p>Read minimum sampling rate and averaging.</p>	<p>Minimum sampling rate for each PPG measurement 1 to 9.  [byte1 to byte9]</p> <p>byteN:  <b>0x01</b> to <b>0x80</b></p>

				(default 0x04)
0x46 for write. AlgoHub	0x07	0x1C [byte1 byte2]	Write maximum DAC offset. byte1: <b>0x00</b> to <b>0x08</b> : PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A)  byte: <b>0x00</b> : 0 uA <b>0x01</b> : 8 uA <b>0x02</b> : 16 uA <b>0x03</b> : 24 uA	-
0x47 for read. AlgoHub	0x07	0x1C	Read maximum DAC offset .	Maximum DAC offset for each PPG measurement 1 to 9. [byte1 to byte9]  ByteN: <b>0x00</b> : 0 uA <b>0x01</b> : 8 uA <b>0x02</b> : 16 uA <b>0x03</b> : 24 uA (default)
0x46 for write. AlgoHub	0x07	0x1D [byte1 MSB LSB]	Write minimum LED current. byte1: <b>0x00</b> to <b>0x08</b> : PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A)  MSB LSB: 16-bit unsigned, lsb=0.1mA	-
0x47 for read. AlgoHub	0x07	0x1D	Read minimum LED current .	Minimum LED current for each PPG measurement 1 to 9. [MSBLSB_1 to MSBLSB_9]  MSBLSB_N: 16-bit unsigned, lsb 0.1mA
0x46 for write. AlgoHub	0x07	0x1E [byte1 MSB LSB]	Write maximum LED current. byte1: <b>0x00</b> to <b>0x08</b> : PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A)  MSB LSB: 16-bit unsigned, lsb=0.1mA	-
0x47 for read. AlgoHub	0x07	0x1E	Read maximum LED current .	Maximum LED current for each PPG measurement 1 to 9. [MSBLSB_1 to MSBLSB_9]  MSBLSB_N: 16-bit unsigned, lsb 0.1mA
0x46 for write. AlgoHub	0x07	0x1F [byte1 MSB LSB]	Write minimum LED current step that the algorithm can increase/decrease with each measurement.	-

			byte1: <b>0x00 to 0x08:</b> PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A)  MSB LSB: 16-bit unsigned, lsb=0.1mA	
0x47 for read. AlgoHub	0x07	0x1F	Read minimum LED current step that the algorithm can increase/decrease with each measurement.	Minimum LED current step for each PPG measurement 1 to 9. [MSBLSB_1 to MSBLSB_9]  MSBLSB_N: 16-bit unsigned, lsb 0.1mA
0x46 for write. AlgoHub	0x07	0x20 [byte1 byte 2]	Write HR master channel selection option for each measurement. byte1: <b>0x00 to 0x08:</b> PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A) byte2: <b>0x00:</b> PD1 selected as master channel. <b>0x01:</b> PD2 selected as master channel. <b>0x02:</b> master channel selected automatically.	-
0x47 for read. AlgoHub	0x07	0x20	Read HR master channel selection option for each measurement .	HR master channel selection for each PPG measurement 1 to 9. [byte1 to byte9]  byteN: <b>0x00:</b> PD1 selected as master channel. <b>0x01:</b> PD2 selected as master channel. <b>0x02:</b> master channel selected automatically.
0x46 for write. AlgoHub	0x07	0x21 [byte1 byte2]	Write full scale PD current. byte1: <b>0x00 to 0x08:</b> PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A) byte2: <b>0x00:</b> 4 uA <b>0x01:</b> 8 uA <b>0x02:</b> 16 uA <b>0x03:</b> 32 uA	-
0x47 for read. AlgoHub	0x07	0x21	Read full scale PD current.	Full scale PD current for each PPG measurement 1 to 9. [byte1 to byte9]  byteN: <b>0x00:</b> 4 uA <b>0x01:</b> 8 uA <b>0x02:</b> 16 uA

				<b>0x03: 32 uA</b>
0x46 for write. AlgoHub	0x07	0x22 [byte]	Write AFE type. byte: <b>0x00</b> : AFE that allows independent integration time and sampling average settings for different measurements, e.g., MAX86171. <b>0x01</b> : AFE that uses shared integration time and sampling average settings for different measurements, e.g., MAX86141	-
0x47 for read. AlgoHub	0x07	0x22	Read AFE type.	AFE type. <b>0x00</b> : AFE that allows independent integration time and sampling average settings for different measurements, e.g., MAX86171. (default) <b>0x01</b> : AFE that uses shared integration time and sampling average settings for different measurements, e.g., MAX86141
0x46 for write. AlgoHub	0x07	0x23 [byte1 byte2]	Write DAC offset for PD channel 1. byte1: <b>0x00</b> to <b>0x08</b> : PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A)  byte2: <b>0x00</b> : 0 uA <b>0x01</b> : 8 uA <b>0x02</b> : 16 uA <b>0x03</b> : 24 uA	-
0x47 for read. AlgoHub	0x07	0x23	Read DAC offset for PD channel 1.	DAC offset for PD channel 1 for each PPG measurement 1 to 9. [byte1 to byte9]  ByteN: <b>0x00</b> : 0 uA (default) <b>0x01</b> : 8 uA <b>0x02</b> : 16 uA <b>0x03</b> : 24 uA
0x46 for write. AlgoHub	0x07	0x24 [byte1 byte2]	Write DAC offset for PD channel 2. byte1: <b>0x00</b> to <b>0x08</b> : PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A)  byte: <b>0x00</b> : 0 uA <b>0x01</b> : 8 uA <b>0x02</b> : 16 uA <b>0x03</b> : 24 uA	-

0x47 for read. AlgoHub	0x07	0x24	Read DAC offset for PD channel 2.	DAC offset for PD channel 2 for each PPG measurement 1 to 9. [byte1 to byte9]  ByteN: <b>0x00</b> : 0 uA (default) <b>0x01</b> : 8 uA <b>0x02</b> : 16 uA <b>0x03</b> : 24 uA
0x46 for write. AlgoHub	0x07	0x25 [byte1 [MSB LSB]	Write initial LED current. (16-bit unsigned, lsb = 0.1mA) byte1: <b>0x00</b> to <b>0x08</b> : PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A)	-
0x47 for read. AlgoHub	0x07	0x25	Read initial LED current.	Initial LED current for each PPG measurement 1 to 9. [MSB1 LSB1 to MSB9 LSB9]  MSBn LSBn: Initial LED current. (16-bit unsigned, lsb = 0.1mA)
0x46 for write. AlgoHub	0x07	0x26 CMD_DELAY = 25ms	Reset AFE and algorithm to the default settings below: 0x50 0x07 0x00: A=-4.4027606, B=18.695252, C= 110.6986082 0x50 0x07 0x0A: 0x00 0x50 0x07 0x06: 0xAF 0x50 0x07 0x07: 0xE4 0x50 0x07 0x13: 0x00 0x50 0x07 0x15: 0x03 0x50 0x07 0x1C: 0x03	-
0x47 for read. AlgoHub	0x07	0x27 (CMD_DELAY = 5 ms)	Read algorithm AFE settings requests.	[byte1 byte2 byte3 byte4 byte5]  byte1, byte2: bit 15: <b>0</b> : No LED current update request <b>1</b> : LED current update requested bits 14-0: LED current (15-bit unsigned, 0.1mA)  byte3: bit 7: <b>0</b> : No integration time update request <b>1</b> : Integration time update requested



				bits 6-0: AFE integration time: <b>0x00</b> : 14.8µs <b>0x01</b> : 29.4µs <b>0x02</b> : 58.7µs <b>0x03</b> : 117.3µs  byte4: bit 7: <b>0</b> : No sampling rate update request <b>1</b> : Sampling rate update requested bits 6-0: Sampling rate and averaging. <b>0x00</b> : 25sps, avg = 1 <b>0x01</b> : 50sps, avg = 2 <b>0x02</b> : 100sps, avg = 4 <b>0x03</b> : 200sps, avg = 8 <b>0x04</b> : 400sps, avg = 16  byte5: DAC offset bit 7: <b>0</b> : No DAC offset update request <b>1</b> : DAC offset update requested <b>0x00</b> : 0 uA <b>0x01</b> : 8 uA <b>0x02</b> : 16 uA <b>0x03</b> : 24 uA
0x46 for write. AlgoHub	0x07	0x28 (CMD_DELAY = 5 ms)	Clear algorithm AFE settings request flag. This command is sent after the host has honored the algorithm AFE requests.	-
0x50 for write. Sensor Hub	0x08	0x00 [A_MSB ... A_LSB] [B_MSB ... B_LSB] [C_MSB ... C_LSB]	Write SpO2 calibration coefficients (12 bytes comprised of three 32-bit signed values, scaled up by 100,000)	-
0x51 for read. Sensor Hub	0x08	0x00	Read SpO2 calibration coefficients (12 bytes comprised of three 32-bit signed values, scaled up by 100,000) Default: 0xFFFF9482C FFE3792B 00A8E9A4 A = -4.4027606 (0xFFFF9482C) B = -18.695252 (0x FFE3792 ) C = 110.6986082 (0x00A8E9A4)	32-bit signed integer A, 32-bit signed integer B, 32-bit signed integer C  Values scaled up by 100,000

0x50 for write. Sensor Hub	0x08	0x01 [MSB LSB]	Write SpO <sub>2</sub> motion-detection period (unsigned 16-bit int, seconds). The algorithm will consider the state to be motionless if the motion is below the threshold for this duration of time.	-
0x51 for read. Sensor Hub	0x08	0x01	Read SpO <sub>2</sub> motion-detection period. The algorithm will consider the state to be motionless if the motion is below the threshold for this duration of time. Default: 0x0002	MSB of period, LSB of period (16-bit unsigned integer, seconds)
0x50 for write. Sensor Hub	0x08	0x02 [MSB .. LSB]	Write SpO <sub>2</sub> motion-detection threshold (signed 32-bit int, equal to 10 <sup>5</sup> x milli-g threshold value)	-
0x51 for read. Sensor Hub	0x08	0x02	Read SpO <sub>2</sub> motion-detection threshold (signed 32-bit int, equal to 10 <sup>5</sup> x milli-g threshold value) Default: 0x01C9C380 (0.3g)	4 bytes (32-bit signed integers which are the milli-g motion threshold times 100,000)
0x50 for write. Sensor Hub	0x08	0x03 [byte]	Write SpO <sub>2</sub> AGC Timeout (sec) .	
0x51 for read. Sensor Hub	0x08	0x03	Read SpO <sub>2</sub> AGC Timeout (sec) . Default: 0x1E	SpO <sub>2</sub> AGC Timeout (8-bit unsigned, seconds)
0x50 for write. Sensor Hub	0x08	0x04	Write the timeout duration for SpO <sub>2</sub> measurement in seconds (1 byte).	-
0x51 for read. Sensor Hub	0x08	0x04	Read the timeout duration for SpO <sub>2</sub> measurement in seconds. Default: 0x78	SpO <sub>2</sub> algorithm timeout (8-bit unsigned, seconds)
0x50 for write. Sensor Hub	0x08	0x05 [byte]	Write initial HR algorithm value (8-bit unsigned).	-
0x51 for read. Sensor Hub	0x08	0x05	Read initial HR algorithm value. Default: 0x3C	Initial heart rate (8-bit unsigned)
0x50 for write. Sensor Hub	0x08	0x06 [MSB] [LSB]	Write height (16-bit unsigned, cm).	-
0x51 for read.	0x08	0x06	Read height. Default: 0x00AF	Height (16-bit unsigned, cm)

Sensor Hub				
0x50 for write. Sensor Hub	0x08	0x07 [MSB] [LSB]	Write weight (16-bit unsigned, kg).	-
0x51 for read. Sensor Hub	0x08	0x07	Read weight. Default: 0x004E	Weight (16-bit unsigned, kg)
0x50 for write. Sensor Hub	0x08	0x08 [byte]	Write age (8-bit unsigned, years).	
0x51 for read. Sensor Hub	0x08	0x08	Read age. Default: 0x1E	Age (8-bit unsigned, years)
0x50 for write Sensor Hub	0x08	0x09 [byte]	Write gender. Byte: <b>0x00</b> : Male <b>0x01</b> : Female	
0x51 for read Sensor Hub	0x08	0x09	Read gender	Gender <b>0x00</b> : Male (default) <b>0x01</b> : Female
0x50 for write. Sensor Hub	0x08	0x0A [WAS algo operation mode byte]	Set the WAS algorithm operation mode (can be switched in runtime): <b>0x00</b> : Continuous HRM + Continuous SpO2 (default). <b>0x01</b> : Continuous HRM + One-Shot SpO2 <b>0x02</b> : Continuous HRM <b>0x03</b> : Sampled HRM <b>0x04</b> : Sampled HRM + One-Shot SpO2 <b>0x05</b> : Activity Tracking ONLY <b>0x06</b> : SpO2 Calibration Data Collection	-
0x51 for read. Sensor Hub	0x08	0x0A	Read the WAS algorithm operation mode.	<b>0x00</b> : Continuous HRM, continuous SpO2 (default) <b>0x01</b> : Continuous HRM, one-shot SpO2 <b>0x02</b> : Continuous HRM <b>0x03</b> : Sampled HRM <b>0x04</b> : Sampled HRM, one-shot SpO2 <b>0x05</b> : Activity tracking only <b>0x06</b> : SpO2 calibration
0x50 for write. Sensor Hub	0x08	0x0B [byte]	Write the enable/disable AEC byte: <b>0x00</b> : Disable <b>0x01</b> : Enable	-

0x51 for read. Sensor Hub	0x08	0x0B	Read the enable/disable AEC.	<b>0x00:</b> Disabled <b>0x01:</b> Enabled (default)
0x50 for write. Sensor Hub	0x08	0x0C [byte]	Write the Skin Contact Detection (SCD) algorithm enable <b>0x00:</b> Disable <b>0x01:</b> Enable	-
0x51 for read. Sensor Hub	0x08	0x0C	Read the Skin Contact Detection (SCD) algorithm enable	<b>0x00:</b> Disable <b>0x01:</b> Enable (default)
0x50 for write. Sensor Hub	0x08	0x0D [MSB] [LSB]	Write adjusted target PD current period (16-bit unsigned, seconds )	-
0x51 for read. Sensor Hub	0x08	0x0D	Read adjusted target PD current period in seconds. Default: 0x0708	Adjusted target PD current period (16-bit unsigned, seconds)
0x50 for write. Sensor Hub	0x08	0x0E [byte]	Write HR motion magnitude threshold (16-bit unsigned, 0.001g)	-
0x51 for read. Sensor Hub	0x08	0x0E [MSB] [LSB]	Read HR motion magnitude threshold. Default: 0x0032 (0.05g)	[MSB] [LSB] motion magnitude threshold (16-bit unsigned, 0.001g)
0x50 for write. Sensor Hub	0x08	0x0F [MSB] [LSB]	Write minimum PD current (16-bit unsigned, 0.1uA).	-
0x51 for read. Sensor Hub	0x08	0x0F	Read minimum PD current. Default: 0x0032	Minimum PD current (16-bit unsigned, 0.1uA)
0x50 for write. Sensor Hub	0x08	0x10 [MSB] [LSB]	Write initial PD current (16-bit unsigned, 0.1uA). This sets the target PD current you would like AEC algorithm to maintain initially. It does not correspond to any register. Once you set what PD current you need, algorithm will calculate the appropriate LED current.	-
0x51 for read. Sensor Hub	0x08	0x10	Read initial PD current. Default: 0x0064	Initial PD current (16-bit unsigned, 0.1uA)
0x50 for write.	0x08	0x11 [MSB] [LSB]	Write target PD current (16-bit unsigned, 0.1uA). Applicable only if Auto Target PD Current Calculation is enabled.	

Sensor Hub				
0x51 for read. Sensor Hub	0x08	0x11	Read target PD current. Default: 0x0064	Target PD current (16-bit unsigned, 0.1uA)
0x50 for write. Sensor Hub	0x08	0x12 [byte]	Write enable/disable automatic calculation of target PD current. byte: <b>0x00</b> : Disable <b>0x01</b> : Enable	-
0x51 for read. Sensor Hub	0x08	0x12	Read enable/disable automatic calculation of target PD current	<b>0x00</b> : Disable <b>0x01</b> : Enable (default)
0x50 for write. Sensor Hub	0x08	0x13 [byte1 byte2]	Write minimum AFE integration time. byte1: <b>0x00</b> to <b>0x08</b> : PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A)  byte2: <b>0x00</b> : 14.8μs <b>0x01</b> : 29.4μs <b>0x02</b> : 58.7μs <b>0x03</b> : 117.3μs	-
0x51 for read. Sensor Hub	0x08	0x13	Read minimum AFE integration time.	Minimum AFE integration time for each PPG measurement 1 to 9. [byte1 to byte9]  byteN: <b>0x00</b> : 14.8μs (default) <b>0x01</b> : 29.4μs <b>0x02</b> : 58.7μs <b>0x03</b> : 117.3μs
0x50 for write. Sensor Hub	0x08	0x14 [byte1 byte2]	Write minimum sampling rate and averaging. byte1: <b>0x00</b> to <b>0x08</b> : PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A)  byte2: <b>0x00</b> : 25sps, avg = 1 <b>0x01</b> : 50sps, avg = 2 <b>0x02</b> : 100sps, avg = 4 <b>0x03</b> : 200sps, avg = 8 <b>0x04</b> : 400sps, avg = 16	-
0x51 for read. Sensor Hub	0x08	0x14	Read minimum sampling rate and averaging.	Minimum sampling rate and averaging for each PPG measurement 1 to 9. [byte1 to byte9]  byteN:

				<b>0x00:</b> 25sps, avg = 1 <b>0x01:</b> 50sps, avg = 2 <b>0x02:</b> 100sps, avg = 4 <b>0x03:</b> 200sps, avg = 8 <b>0x04:</b> 400sps, avg = 16
0x50 for write. Sensor Hub	0x08	0x15 [byte1 byte2]	Write maximum integration time: byte1: <b>0x00 to 0x08:</b> PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A) byte2: <b>0x00:</b> 14.8μs <b>0x01:</b> 29.4μs <b>0x02:</b> 58.7μs <b>0x03:</b> 117.3μs	-
0x51 for read. Sensor Hub	0x08	0x15	Read maximum integration time:	Maximum integration time for each PPG measurement 1 to 9. [byte1 to byte9] byteN: <b>0x00:</b> 14.8μs <b>0x01:</b> 29.4μs <b>0x02:</b> 58.7μs <b>0x03:</b> 117.3μs (default)
0x50 for write. Sensor Hub	0x08	0x16 [byte1 byte2]	Write maximum sampling rate and averaging: byte1: <b>0x00 to 0x08:</b> PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A) byte2: <b>0x00:</b> 25sps, avg = 1 <b>0x01:</b> 50sps, avg = 2 <b>0x02:</b> 100sps, avg = 4 <b>0x03:</b> 200sps, avg = 8 <b>0x04:</b> 400sps, avg = 16	-
0x51 for read. Sensor Hub	0x08	0x16	Read maximum sampling rate and averaging:	Maximum sampling rate and averaging for each PPG measurement 1 to 9. [byte1 to byte9] byteN: <b>0x00:</b> 25sps, avg = 1 <b>0x01:</b> 50sps, avg = 2 <b>0x02:</b> 100sps, avg = 4 <b>0x03:</b> 200sps, avg = 8 <b>0x04:</b> 400sps, avg = 16 (default)
0x50 for write. Sensor Hub	0x08	0x17 0xWX, 0xYZ	Write slot and PD configuration for the two HR inputs to the WHRM algorithm. 0xWX, 0xYZ:	-



			<p>WX is input 1 of the WHRM algorithm. W (MS nibble of WX): W = 0 for Slot 1 W = 1 for Slot 2 W = 2 for Slot 3 W = 3 for Slot 4 W = 4 for Slot 5 W = 5 for Slot 6 W = 7 for Slot not used X (LS nibble of WX): X = 0 for PD1 X = 1 for PD2 X = 3 for PD not used.</p> <p>YZ is input 2 of the WHRM algorithm. Y (MS nibble of YZ): Y = 0 for Slot 1 Y = 1 for Slot 2 Y = 2 for Slot 3 Y = 3 for Slot 4 Y = 4 for Slot 5 Y = 5 for Slot 6 Y = 7 for Slot not used Z (LS nibble of YZ): Z = 0 for PD1 Z = 1 for PD2 Z = 3 for PD not used.</p>	
0x51 for read. Sensor Hub	0x08	0x17	<p>Read Slot and PD configuration for the two HR inputs to the WHRM algorithm.</p>	<p>0xWX, 0xYZ Default: 0x0001</p> <p>WX is input 1 of the WHRM algorithm. W (MS nibble of WX): W = 0 for Slot 1 W = 1 for Slot 2 W = 2 for Slot 3 W = 3 for Slot 4 W = 4 for Slot 5 W = 5 for Slot 6 W = 7 for Slot not used X (LS nibble of WX): X = 0 for PD1 X = 1 for PD2 X = 3 for PD not used.</p> <p>YZ is input 2 of the WHRM algorithm. Y (MS nibble of YZ): Y = 0 for Slot 1 Y = 1 for Slot 2 Y = 2 for Slot 3 Y = 3 for Slot 4 Y = 4 for Slot 5 Y = 5 for Slot 6</p>

				Y = 7 for Slot not used Z (LS nibble of YZ): Z = 0 for PD1 Z = 1 for PD2 Z = 3 for PD not used.
0x50 for write. Sensor Hub	0x08	0x18 0xWX 0xYZ	Write Slot and PD configuration for the IR, red inputs to the WSpO2 algorithm.  0xWX, 0xYZ:  WX is the LED/PD used for IR for the WSpO2 algorithm. W (MS nibble of WX): W = 0 for Slot 1 W = 1 for Slot 2 W = 2 for Slot 3 W = 3 for Slot 4 W = 4 for Slot 5 W = 5 for Slot 6 W = 7 for Slot not used X (LS nibble of WX): X = 0 for PD1 X = 1 for PD2 X = 3 for PD not used.  YZ is the LED/PD used for red for the WSpO2 algorithm. Y (MS nibble of YZ): Y = 0 for Slot 1 Y = 1 for Slot 2 Y = 2 for Slot 3 Y = 3 for Slot 4 Y = 4 for Slot 5 Y = 5 for Slot 6 Y = 7 for Slot not used Z (LS nibble of YZ): Z = 0 for PD1 Z = 1 for PD2 Z = 3 for PD not used	-
0x51 for read. Sensor Hub	0x08	0x18	Read Slot and PD configuration for the IR, red inputs to the WSpO2 algorithm.	0xWX, 0xYZ Default: 0x1020  WX is the LED/PD used for IR for the WSpO2 algorithm. W (MS nibble of WX): W = 0 for Slot 1 W = 1 for Slot 2 W = 2 for Slot 3 W = 3 for Slot 4 W = 4 for Slot 5 W = 5 for Slot 6 W = 7 for Slot not used X (LS nibble of WX):

				<p>X = 0 for PD1 X = 1 for PD2 X = 3 for PD not used.</p> <p>YZ is the LED/PD used for red for the WSpO<sub>2</sub> algorithm. Y (MS nibble of YZ): Y = 0 for Slot 1 Y = 1 for Slot 2 Y = 2 for Slot 3 Y = 3 for Slot 4 Y = 4 for Slot 5 Y = 5 for Slot 6 Y = 7 for Slot not used Z (LS nibble of YZ): Z = 0 for PD1 Z = 1 for PD2 Z = 3 for PD not used</p>
0x50 for write. Sensor Hub	0x08	0x19 [byte1 WX]	<p>Write LED driver and MUX configuration. byte1: <b>0x00 to 0x08</b>: PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A)</p> <p>W (MS nibble of WX): W = 0 for LED Driver B will fire LED1 W = 1 for LED Driver B will fire LED2 W = 2 for LED Driver B will fire LED3 W = 3 for LED Driver B will fire LED4 W = 4 for LED Driver B will fire LED5 W = 5 for LED Driver B will fire LED6</p> <p>X (LS nibble of WX): X = 0 for LED Driver A will fire LED1 X = 1 for LED Driver A will fire LED2 X = 2 for LED Driver A will fire LED3 X = 3 for LED Driver A will fire LED4 X = 4 for LED Driver A will fire LED5 X = 5 for LED Driver A will fire LED6</p>	
0x51 for read. Sensor Hub	0x08	0x19	<p>Read LED driver and MUX .</p>	<p>LED driver and MUX configuration for each PPG measurement 1 to 9. [WX1 to WX9]</p> <p>W (MS nibble of WX_N): W = 0 for LED Driver B will fire LED1 W = 1 for LED Driver B will fire LED2 W = 2 for LED Driver B will fire LED3 W = 3 for LED Driver B will fire LED4</p>

				W = 4 for LED Driver B will fire LED5 W = 5 for LED Driver B will fire LED6  X (LS nibble of WX_N): X = 0 for LED Driver A will fire LED1 X = 1 for LED Driver A will fire LED2 X = 2 for LED Driver A will fire LED3 X = 3 for LED Driver A will fire LED4 X = 4 for LED Driver A will fire LED5 X = 5 for LED Driver A will fire LED6
0x50 for write. Sensor Hub	0x08	0x1A [byte1 byte2]	Write initial AFE integration time. byte1: <b>0x00</b> to <b>0x08</b> : PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A)  byte2: <b>0x00</b> : 14.8μs <b>0x01</b> : 29.4μs <b>0x02</b> : 58.7μs <b>0x03</b> : 117.3μs	-
0x51 for read. Sensor Hub	0x08	0x1A	Read initial AFE integration time.	Initial AFE integration time for each PPG measurement 1 to 9. [byte1 to byte9]  ByteN: <b>0x00</b> : 14.8μs <b>0x01</b> : 29.4μs <b>0x02</b> : 58.7μs <b>0x03</b> : 117.3μs (default)
0x50 for write. Sensor Hub	0x08	0x1B [byte1 byte2]	Write initial sampling rate and averaging. <b>0x00</b> to <b>0x08</b> : PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A)  byte2: <b>0x00</b> : 25sps, avg = 1 <b>0x01</b> : 50sps, avg = 2 <b>0x02</b> : 100sps, avg = 4 <b>0x03</b> : 200sps, avg = 8 <b>0x04</b> : 400sps, avg = 16	
0x51 for read. Sensor Hub	0x08	0x1B	Read minimum sampling rate and averaging.	Minimum sampling rate and averaging for each PPG measurement 1 to 9. [byte1 to byte9]  ByteN:

				<b>0x00:</b> 25sps, avg = 1 <b>0x01:</b> 50sps, avg = 2 <b>0x02:</b> 100sps, avg = 4 (default) <b>0x03:</b> 200sps, avg = 8 <b>0x04:</b> 400sps, avg = 16
0x50 for write. Sensor Hub	0x08	0x1C [byte1 byte2]	Write maximum DAC offset. byte1: <b>0x00 to 0x08:</b> PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A)  byte: <b>0x00:</b> 0 uA <b>0x01:</b> 8 uA <b>0x02:</b> 16 uA <b>0x03:</b> 24 uA	-
0x51 for read. Sensor Hub	0x08	0x1C	Read maximum DAC offset .	Maximum DAC offset for each PPG measurement 1 to 9. [byte1 to byte9]  ByteN: <b>0x00:</b> 0 uA <b>0x01:</b> 8 uA <b>0x02:</b> 16 uA <b>0x03:</b> 24 uA (default)
0x50 for write. Sensor Hub	0x08	0x1D [byte1 MSB LSB]	Write minimum LED current. byte1: <b>0x00 to 0x08:</b> PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A)  MSB LSB: 16-bit unsigned, lsb=0.1mA	-
0x51 for read. Sensor Hub	0x08	0x1D	Read minimum LED current .	Minimum LED current for each PPG measurement 1 to 9. [MSBLSB_1 to MSBLSB_9]  MSBLSB_N: 16-bit unsigned, lsb 0.1mA
0x50 for write. Sensor Hub	0x08	0x1E [byte1 MSB LSB]	Write maximum LED current. byte1: <b>0x00 to 0x08:</b> PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A)  MSB LSB: 16-bit unsigned, lsb=0.1mA	-
0x51 for read. Sensor Hub	0x08	0x1E	Read maximum LED current .	Maximum LED current for each PPG measurement 1 to 9. [MSBLSB_1 to MSBLSB_9]

				MSBLSB_N: 16-bit unsigned, lsb 0.1mA
0x50 for write. Sensor Hub	0x08	0x1F [byte1 MSB LSB]	Write minimum LED current step that the algorithm can increase/decrease with each measurement. byte1: <b>0x00</b> to <b>0x08</b> : PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A)  MSB LSB: 16-bit unsigned, lsb=0.1mA	-
0x51 for read. Sensor Hub	0x08	0x1F	Read minimum LED current step that the algorithm can increase/decrease with each measurement.	Minimum LED current step for each PPG measurement 1 to 9. [MSBLSB_1 to MSBLSB_9]  MSBLSB_N: 16-bit unsigned, lsb 0.1mA
0x50 for write. Sensor Hub	0x08	0x20 [byte1 byte 2]	Write HR master channel selection option for each measurement. byte1: <b>0x00</b> to <b>0x08</b> : PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A) byte2: <b>0x00</b> : PD1 selected as master channel. <b>0x01</b> : PD2 selected as master channel. <b>0x02</b> : master channel selected automatically.	-
0x51 for read. Sensor Hub	0x08	0x20	Read HR master channel selection option for each measurement .	HR master channel selection for each PPG measurement 1 to 9. [byte1 to byte9]  byteN: <b>0x00</b> : PD1 selected as master channel. <b>0x01</b> : PD2 selected as master channel. <b>0x02</b> : master channel selected automatically.
0x50 for write. Sensor Hub	0x08	0x21 [byte1 byte2]	Write full scale PD current. byte1: <b>0x00</b> to <b>0x08</b> : PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A) byte2: <b>0x00</b> : 4 uA <b>0x01</b> : 8 uA <b>0x02</b> : 16 uA <b>0x03</b> : 32 uA	-
0x51 for read.	0x08	0x21	Read full scale PD current.	Full scale PD current for each PPG measurement 1 to 9.

Sensor Hub				[byte1 to byte9]  byteN: <b>0x00:</b> 4 uA <b>0x01:</b> 8 uA <b>0x02:</b> 16 uA <b>0x03:</b> 32 uA
0x50 for write. Sensor Hub	0x08	0x23 [byte1 byte2]	Write DAC offset for PD channel 1. byte1: <b>0x00 to 0x08:</b> PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A)  byte2: <b>0x00:</b> 0 uA <b>0x01:</b> 8 uA <b>0x02:</b> 16 uA <b>0x03:</b> 24 uA	-
0x51 for read. Sensor Hub	0x08	0x23	Read DAC offset for PD channel 1.	DAC offset for PD channel 1 for each PPG measurement 1 to 9. [byte1 to byte9]  ByteN: <b>0x00:</b> 0 uA (default) <b>0x01:</b> 8 uA <b>0x02:</b> 16 uA <b>0x03:</b> 24 uA
0x50 for write. Sensor Hub	0x08	0x24 [byte1 byte2]	Write DAC offset for PD channel 2. byte1: <b>0x00 to 0x08:</b> PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A)  byte: <b>0x00:</b> 0 uA <b>0x01:</b> 8 uA <b>0x02:</b> 16 uA <b>0x03:</b> 24 uA	-
0x51 for read. Sensor Hub	0x08	0x24	Read DAC offset for PD channel 2.	DAC offset for PD channel 2 for each PPG measurement 1 to 9. [byte1 to byte9]  ByteN: <b>0x00:</b> 0 uA (default) <b>0x01:</b> 8 uA <b>0x02:</b> 16 uA <b>0x03:</b> 24 uA
0x50 for write. Sensor Hub	0x08	0x25 [byte1 [MSB LSB]	Write initial LED current. (16-bit unsigned, lsb = 0.1mA)  byte1: <b>0x00 to 0x08:</b> PPG measurement 1 to 9. (1: Green, 2: IR, 3: red, N/A, ... N/A)	-



0x51 for read. Sensor Hub	0x08	0x25	Read initial LED current.	Initial LED current for each PPG measurement 1 to 9. [MSB1 LSB1 to MSB9 LSB9]  MSBn LSBn: Initial LED current. (16-bit unsigned, lsb = 0.1mA)
0x50 for write. Sensor Hub	0x08	0x40 [byte1]	Write biometric operation mode.  byte1: <b>0x00</b> : Raw <b>0x01</b> : WAS (WHRM, WSpO2) <b>0x02</b> : BPT (Blood Pressure Trending) <b>0x03</b> : BPT+WAS extended (WHRM, WSpO2)	
0x51 for read. Sensor Hub	0x08	0x40	Read biometric operation mode	<b>0x00</b> : Raw <b>0x01</b> : WAS (WHRM, WSpO2) (default) <b>0x02</b> : BPT (Blood Pressure Trending) <b>0x03</b> : BPT+WAS extended (WHRM, WSpO2)
0x52 for write. Sensor Hub	0x08	Definition is dependent on biometric operation mode.  Biometric operation mode = Raw: <b>0x00</b> : Disable (CMD_DELAY = 200ms) <b>0x01</b> : Enable raw report (CMD_DELAY = 500ms)  Biometric operation mode = WAS: <b>0x00</b> : Disable (CMD_DELAY = 200ms) <b>0x01</b> : Enable normal algorithm report (CMD_DELAY = 500ms) <b>0x02</b> : N/A	Write enable/disable the biometric operation mode algorithm, SensorHub.	

		<p>Biometric operation mode = <i>BPT</i>:  <b>0x00</b>: Disable (CMD_DELAY = 200ms)  <b>0x01</b>: Enable BPT  Calibration report (CMD_DELAY = 500ms)  <b>0x02</b>: Enable BPT  Estimation algorithm report (CMD_DELAY = 500ms)</p> <p>Biometric operation mode = <i>BPT+WAS</i>:  <b>0x00</b>: Disable (CMD_DELAY = 200ms)  <b>0x01</b>: Enable BPT+WAS algorithm report (CMD_DELAY = 500ms)  <b>0x02</b>: Enable <i>BPT+WAS extended</i> algorithm report (CMD_DELAY = 500ms)</p>		
0x50 for write. Sensor Hub	0x08	0x51 [byte1]	Set BPT run mode to estimation/calibration [byte1] <b>0x00</b> : BPT Calibration <b>0x01</b> : BPT Estimation	
0x51 for read. Sensor Hub	0x08	0x51 [byte1]	Read the BPT run mode.	<b>0x00</b> : BPT Calibration <b>0x01</b> : BPT Estimation

0x50 for write. Sensor Hub	0x08	0x60, [byte1...byte240]	<p>Write the BPT user calibration vector. To restore user calibration data after a user has restarted the host app:</p> <ul style="list-style-type: none"> <li>• write the cal_index</li> <li>• write the 240 bytes of BPT calibration vector and stored in host memory.</li> </ul> <p>Repeat the bulleted items until max(cal_index)</p>	
0x51 for read. Sensor Hub	0x08	0x60	<p>Read BPT user calibration vector, after calibration is complete (BPTstatus is 2 and BPT progress is 100%) and store in host memory as an array of 240 bytes vectors indexed by cal_index.</p>	<p>[byte1...byte4]: LSB...MSB BPT date</p> <p>[byte5...byte8]: LSB...MSB BPT time</p> <p>[byte9...byte12]: LSB...MSB reference systolic</p> <p>[byte13...byte16]: LSB...MSB reference systolic</p> <p>[byte17...byte20]: LSB...MSB reference systolic</p> <p>[byte21...240]: user calibration data</p>
0x50 for write. Sensor Hub	0x08	0x61	<p>Write the BPTcalibration date and time as two 32-bit unsigned values (little-endian format) in the following order:</p> <p>value1: BPT date in YYYYMMDD decimal form</p> <p>value2: BPT time in HHMMSS decimal form</p>	<p>[value1 .. value2]</p> <p>Example: 0x61, 0x3D, 0x34, 0x01, 0xE0, 0xDF, 0x01, 0x00 for date = 20200801, time = 12:28:48</p>
0x50 for write. Sensor Hub	0x08	0x62 [byte1 byte2 byte3]	<p>Write the cal_index, reference systolic, reference diastolic. This command is sent after the command 50, 08, 04, 02 (BPT operation mode) and before the 05 08 01 (enable algorithm BPT calibration) as part of the multi-point calibration procedure.</p> <p>byte1: <i>cal_index</i>: 0 to 5: Set this to 0 for the first subject calibration.</p> <p>byte2: <i>reference systolic</i></p> <p>byte3: <i>reference diastolic</i></p>	
0x50 for write. Sensor Hub	0x08	0x63 [byte1]	<p>Write the cal_index. Set the cal_index 0 to 5. This command is sent before the command 50 04 03 (set calibration vector command) when restoring the user calibration vectors.</p> <p>byte1: <i>cal_index</i></p>	
0x51 for read.	0x08	0x63	Read the cal_index: byte1.	[byte1]

Sensor Hub				
0x50 for write. Sensor Hub	0x08	0x64 [byte1]	Write the BPT continuous/one-shot mode.  byte1: <b>0x00</b> : One-shot BPT <b>0x01</b> : Continuous BPT	
0x51 for read. Sensor Hub	0x08	0x64	Read the BPT continuous/one-shot mode.	[byte1]  <b>0x00</b> : One-shot BPT <b>0x01</b> : Continuous BPT
0x54 for write	[byte]	-	byte: <b>0x00</b> : MAX32674C releases sensor SPI bus to Host. AlgoHub mode. <b>0x01</b> : sensor SPI bus controlled by MAX32674C. SensorHub mode.	-
0x80 Bootloader	0x00	Use bytes 0x28 to 0x32 from the .msbl file as the IV bytes.	Bootloader mode flash the application .msbl: Set the initialization vector (IV) bytes.	-
0x80 Bootloader	0x01	Use bytes 0x34 to 0x43 from the .msbl file.	Bootloader mode flash the application .msbl: Set the authentication bytes.	-
0x80 Bootloader	0x02	0x00, Number of pages located at byte 0x44 from the .msbl file.	Bootloader mode flash the application .msbl: Set the number of pages to flash.	-
0x80 Bootloader	0x03	- (CMD_DELAY = 1400ms)	Bootloader mode flash the application .msbl: Erase the application flash memory.	-
0x80 Bootloader	0x04	The first page is specified by byte 0x4C from the .msbl file. The total bytes for each message protocol are the page size plus 16 bytes of CRC. (CMD_DELAY = 680ms)	Bootloader mode flash the application .msbl: Send the page values. Each page sent includes 16 CRC bytes for that page, so there are 8208 bytes per page sent in the payload of the message.	-
0x80 Bootloader	0x06	[MSB] [LSB]	Bootloader mode flash the application .msbl: Size of the partial page if	-

			using partial page writing (16-bit unsigned, range 1 to 8208)	
0x81 Bootloader	0x00	-	Bootloader mode: Get bootloader version.	Major version byte, Minor version byte, Revision byte
0x81 Bootloader	0x01	-	Bootloader mode flash the application .msbl: Get the page size in bytes.	Upper byte of page size, Lower byte of page size
0x81 Bootloader	0x02	-	Get unique serial number (USN).	USN: 24 bytes
0x82 Bootloader Configuration	0x00 (CMD_DELAY = 300ms)		Save bootloader configurations. Write this command after changes are made to any of the Bootloader Configuration settings. The device should be restarted for the new configuration to be active.	
0x82 Bootloader Configuration	0x01	0x00	0x01 Configure the device to check the state of the MFIO GPIO pin to decide if entering bootloader mode or application mode on after a reset.	
0x82 Bootloader Configuration	0x01	0x01	x00 Configure the MFIO pin to be P0.0	
0x82 Bootloader Configuration	0x01	0x20	0x00 Configure the MFIO pin the be use active-low for bootloader mode. The device enters bootloader mode if the MFIO GPIO pin is held low during power-on or during a RSTN device pin cycle.	
0xFF	0x03		Read the AlgoHubSensorHub version.	Major version byte, Minor version byte, Revision byte

## Configuring the AlgoHub and Reading Processed Data

### AlgoHub Output FIFO Format

The algorithm output FIFO format for AlgoHub is listed below (The initial response byte, “Read Status Byte” is not shown) . The algorithm outputs 24 bytes to the output FIFO.

**Table 3. AlgoHub Output FIFO Format**

DATA ITEM	# OF BYTES (MSB FIRST)	DESCRIPTION
PPG	PPG1	0,0,0
AFE Settings Request Flag, Algorithm operation mode	1	bit 8: 0: no AFE settings request flag 1: AFE settings request flag set; Host to issue the command 0x47 0x07 0x27 to read the AFE requested settings and then issue 0x47 0x07 0x28 when the request has been honored.  bits 6 to 0: Algorithm operation mode: 0: Continuous HRM and Continuous SpO <sub>2</sub> 1: Continuous HRM and One-Shot SpO <sub>2</sub> 2: Continuous HRM 3: Sampled HRM 4: Sampled HRM and One-Shot SpO <sub>2</sub> 5: Activity tracking 6: SpO <sub>2</sub> Calibration Data Collection
HR	2	10x Calculated heart rate
HR confidence	1	Calculated confidence level in %
RR	2	10x RR – inter-beat interval in ms Only shows a nonzero value when a new value is calculated.
RR confidence	1	Calculated confidence level of RR in % Only shows a nonzero value when a new value is calculated.
Activity class	1	Activity class: 0: Rest 1: Other 2: Walk 3: Run 4: Bike
R	2	1000x Calculated SpO <sub>2</sub> R value
SpO <sub>2</sub> confidence	1	Calculated SpO <sub>2</sub> confidence level in %
SpO <sub>2</sub>	2	10x Calculated SpO <sub>2</sub> %
SpO <sub>2</sub> percent complete	1	Calculation progress in percent in one-shot mode of algorithm. In continuous mode, it is reported as zero and only jumps to 100 when the SpO <sub>2</sub> value is updated.
SpO <sub>2</sub> low signal quality flag	1	Shows the low quality of the PPG signal: 0: Good quality 1: Low quality
SpO <sub>2</sub> motion6001200 flag	1	Shows excessive motion: 0: No motion 1: Excessive motion
SpO <sub>2</sub> low PI flag	1	Shows the low perfusion index (PI) of the PPG signal: 0: Normal PI

		1: Low PI
SpO <sub>2</sub> unreliable R flag	1	Shows the reliability of R: 0: Reliable 1: Unreliable
SpO <sub>2</sub> state	1	Reported status of the SpO <sub>2</sub> algorithm: 0: LED adjustment 1: Computation 2: Success 3: Timeout
SCD state	1	Skin contact state: 0: Undetected 1: Off skin 2: On some subject 3: On skin
Algorithm status	1	Algorithm return status: 0: Success 1: Null pointer error 2: Instance already initialized 3: Uninitialized instance 4: AFE controller error 5: Error, not categorized 6: Incompatible internal library error

### AlgoHub AEC Normal Algorithm Report Configuration and Batched Processing

Automatic Exposure Control (AEC) is Maxim's gain control algorithm that is superior to AGC. The AEC algorithm optimally maintains the best SNR range and power optimization. The targeted SNR range is maintained regardless of skin color or ambient temperature within the limits of the LED currents configurations; The AEC dynamically manages the appropriate register settings for sampling rate, LED current, pulse width and integration time.

For AlgoHub mode, the host is responsible for initializing the AFE, accel registers and informing the AlgoHub of the AFE initialization settings and for updating the AFE registers when the AFE request flag is set.

In the example below, both the AEC and SCD are enabled. The sequence of commands is shown in the table below.



**Table 4. Host Commands—AlgoHub AEC and Read Normal Algorithm Report**

	#	HOST COMMAND (HEX)	COMMAND DESCRIPTION	RESPONSE (HEX)
START ALGORITHM		Host initializes the AlgoHub and starts the algorithm using following commands:		
	1.1	AA 02 00 (optional)	Read the operating mode	AB 00 00 application mode
	1.2	AA FF 03 (optional)	Read the AlgoHubSensorHub version for MAX32674C (x.y.z)	AB 00 32 YY ZZ
	1.3	AA 54 00	MAX32674C releases sensor SPI bus to Host. AlgoHub mode.	AB 00
	1.4	AA 10 01 01	Set the AlgoHub interrupt threshold.	AB 00
	1.5	AA 46 07 0B 01	Enable AEC	AB 00
	1.6	AA 46 07 12 01	Enable automatic PD calculation	AB 00
	1.7	AA 46 07 1A 00 03	Set PPG measurement 1 AFE integration time to 117.3us	AB 00
	1.8	AA 46 07 1A 01 03	Set PPG measurement 2 AFE integration time to 117.3us	AB 00
	1.9	AA 46 07 1A 02 03	Set PPG measurement 3 AFE integration time to 117.3us	AB 00
	1.10	AA 46 07 1B 00 01	Set PPG measurement 1 AFE sampling rate to 50sps, averaging to 2	AB 00
	1.11	AA 46 07 1B 01 04	Set PPG measurement 2 AFE sampling rate to 400sps, averaging to 16	AB 00
	1.12	AA 46 07 1B 02 04	Set PPG measurement 3 AFE sampling rate to 400sps, averaging to 16	AB 00
	1.13	AA 46 07 0F 00 7D	Set minimum PD current to 12.5uA	AB 00
	1.14	AA 46 07 10 01 38	Set initial PD current to 31.2uA	AB 00
	1.15	AA 46 07 0D 07 08	Set the adjusted target PD current period to 1800s	AB 00
	1.16	AA 46 07 0E 00 32	Set the HR motion magnitude threshold to 0.05g.	AB 00
	1.17	AA 46 07 23 00 00	Set the DAC offset for PD channel 1 for PPG measurement 1	AB 00
	1.18	AA 46 07 24 00 00	Set the DAC offset for PD channel 2 for PPG measurement 1	AB 00
	1.19	AA 46 07 1C 00 00	Set the maximum DAC offset for PPG measurement 1 to 0	AB 00
	1.20	AA 46 07 23 01 00	Set the DAC offset for PD channel 1 for PPG measurement 2	AB 00
	1.21	AA 46 07 24 01 00	Set the DAC offset for PD channel 2 for PPG measurement 2	AB 00
	1.22	AA 46 07 23 02 00	Set the DAC offset for PD channel 3 for PPG measurement 1	AB 00
	1.23	AA 46 07 24 02 00	Set the DAC offset for PD channel 2 for PPG measurement 3	AB 00
	1.24	AA 46 07 11 01 38	Set target PD current to 31.2uA	AB 00
	1.25	AA 46 07 0C 01	Enable SCD	AB 00
	1.26	AA 46 07 25 00 00 64	Set initial LED current to 10mA for PPG measurement 1	AB 00
	1.27	AA 46 07 25 01 00 C8	Set initial LED current to 20mA for PPG measurement 2	AB 00

READING SAMPLES REPORT IN OUTPUT FIFO	1.28	AA 46 07 25 02 00 C8	Set initial LED current to 20mA for PPG measurement 3	AB 00
	1.29	AA 10 00 03	Set the output of the AlgoHub to sensor and algorithm.	AB 00
	1.30	AA 44 07 01 01 (CMD_DELAY = 500ms)	Enable Wearable Algorithm Suite AlgoHub	AB 00
	Host repeats the following commands every 40ms			
	2.1	AA 14 00 [PPG1 PPG2 PPG3 PPG4 PPG5 PPG6 ACCLX ACCLY ACCLZ] (CMD_DELAY = 16 ms)	Write LEDPPG, accel data to the algorithm input FIFO. PPG1: 3 bytes [MSB .. LSB] WHRM HR Channel 1 (green1) PPG2: 3 bytes [MSB .. LSB] WHRM HR Channel 2 (green2) PPG3: 3 bytes [MSB .. LSB] SpO2 IR Channel (IR) PPG4: 3 bytes [MSB .. LSB] SpO2 IR Channel (red) PPG5: 3 bytes [MSB .. LSB] PPG6: 3 bytes [MSB .. LSB] Accelerometer set to $\pm 8g$ : ACCLX: 2 byte accel X value (2's complement .001g): [MSB LSB] ACCLY: 2 byte accel Y value (2's complement .001g): [MSB LSB] ACCLZ: 2 byte accel Z value (2's complement .001g): [MSB LSB]	AB 00 mm nn  mmnn is the ppg, accel number of bytes received (16-bit unsigned)
		Delay for 20 msec		
	2.2	AA 00 00	Read the AlgoHub status byte: Bit 0: Sensor comm error Bits 1 and 2: Reserved Bit 3: FIFO filled to threshold (DataRdyInt) Bit 4: Output FIFO overflow (FifoOutOvrInt) Bit 5: Input FIFO overflow (FifoInOvrInt) Bit 6: AlgoHub busy (DevBusy) Bit 7: Reserved If DataRdyInt is set, proceed to the next step.	AB 00 08
	2.3	AA 12 01 (CMD_DELAY = 5 ms)	Read the data stored in the FIFO. The format of the samples report is shown in the normal algorithm report table.	AB 00 [Normal Algorithm Report]
	2.4, 2.5 are commanded if the AFE settings request flag is set			
	2.4	AA 47 07 27 (CMD_DELAY = 5 ms)	Read algorithm AFE settings requests.	[byte1 byte2 byte3 byte4 byte5]  byte1, byte2: bit 15: 0: No LED current update request 1: LED current update requested bits 14-0: LED current (15-bit unsigned, 0.1mA)  byte3: bit 7:

				<p><b>0:</b> No integration time update request  <b>1:</b> Integration time update requested  bits 6-0:  AFE integration time:  <b>0x00:</b> 14.8µs  <b>0x01:</b> 29.4µs  <b>0x02:</b> 58.7µs  <b>0x03:</b> 117.3µs</p> <p>byte4:  bit 7:  <b>0:</b> No sampling rate update request  <b>1:</b> Sampling rate update requested  bits 6-0:  Sampling rate and averaging.  <b>0x00:</b> 25sps, avg = 1  <b>0x01:</b> 50sps, avg = 2  <b>0x02:</b> 100sps, avg = 4  <b>0x03:</b> 200sps, avg = 8  <b>0x04:</b> 400sps, avg = 16</p> <p>byte5: DAC offset  bit 7:  <b>0:</b> No DAC offset update request  <b>1:</b> DAC offset update requested  <b>0x00:</b> 0 µA  <b>0x01:</b> 8 µA  <b>0x02:</b> 16 µA  <b>0x03:</b> 24 µA</p>
	2.5	AA 47 07 28 (CMD_DELAY = 5 ms)	Clear algorithm AFE settings request flag. This command is sent after the host has honored the algorithm AFE requests.	AB 00 08
STOP	Host disables sensors and algorithm:			
	3.1	AA 44 07 00 01 (CMD_DELAY = 200ms)	Disable the WAS algorithm.	AB 00
	3.2	AA 46 26 (CMD_DELAY = 25ms)	Reset the AFE settings of the WAS algorithm.	AB 00

The example below illustrates how the host may utilize batch processing to communicate with the AlgoHub at a 1s rate. The sequence of commands is shown in the table below.

**Table 5. Host Commands—Batched Processing AlgoHub AEC and Read Normal Algorithm Report**

	#	HOST COMMAND (HEX)	COMMAND DESCRIPTION	RESPONSE (HEX)
START ALGORITHM		Host initializes the AlgoHub and starts the algorithm using following commands:		
	1.1	AA 02 00 (optional)	Read the operating mode	AB 00 00 application mode
	1.2	AA FF 03 (optional)	Read the AlgoHubSensorHub version for MAX32674C (x.y.z)	AB 00 32 YY ZZ
	1.3	AA 54 00	MAX32674C releases sensor SPI bus to Host. AlgoHub mode.	AB 00
	1.4	AA 10 01 01	Set the AlgoHub interrupt threshold.	AB 00
	1.5	AA 46 07 0B 01	Enable AEC	AB 00
	1.6	AA 46 07 12 01	Enable automatic PD calculation	AB 00
	1.7	AA 46 07 1A 00 03	Set PPG measurement 1 AFE integration time to 117.3us	AB 00
	1.8	AA 46 07 1A 01 03	Set PPG measurement 2 AFE integration time to 117.3us	AB 00
	1.9	AA 46 07 1A 02 03	Set PPG measurement 3 AFE integration time to 117.3us	AB 00
	1.10	AA 46 07 1B 00 01	Set PPG measurement 1 AFE sampling rate to 50sps, averaging to 2	AB 00
	1.11	AA 46 07 1B 01 04	Set PPG measurement 2 AFE sampling rate to 400sps, averaging to 16	AB 00
	1.12	AA 46 07 1B 02 04	Set PPG measurement 3 AFE sampling rate to 400sps, averaging to 16	AB 00
	1.13	AA 46 07 0F 00 7D	Set minimum PD current to 12.5uA	AB 00
	1.14	AA 46 07 10 01 38	Set initial PD current to 31.2uA	AB 00
	1.15	AA 46 07 0D 07 08	Set the adjusted target PD current period to 1800s	AB 00
	1.16	AA 46 07 0E 00 32	Set the HR motion magnitude threshold to 0.05g.	AB 00
	1.17	AA 46 07 23 00 00	Set the DAC offset for PD channel 1 for PPG measurement 1	AB 00
	1.18	AA 46 07 24 00 00	Set the DAC offset for PD channel 2 for PPG measurement 1	AB 00
	1.19	AA 46 07 1C 00 00	Set the maximum DAC offset for PPG measurement 1 to 0	AB 00
	1.20	AA 46 07 23 01 00	Set the DAC offset for PD channel 1 for PPG measurement 2	AB 00
	1.21	AA 46 07 24 01 00	Set the DAC offset for PD channel 2 for PPG measurement 2	AB 00
	1.22	AA 46 07 23 02 00	Set the DAC offset for PD channel 3 for PPG measurement 1	AB 00
	1.23	AA 46 07 24 02 00	Set the DAC offset for PD channel 2 for PPG measurement 3	AB 00
	1.24	AA 46 07 11 01 38	Set target PD current to 31.2uA	AB 00
	1.25	AA 46 07 0C 01	Enable SCD	AB 00
	1.26	AA 46 07 25 00 00 64	Set initial LED current to 10mA for PPG measurement 1	AB 00
	1.27	AA 46 07 25 01 00 C8	Set initial LED current to 20mA for PPG measurement 2	AB 00
	1.28	AA 46 07 25 02 00 C8	Set initial LED current to 20mA for PPG measurement 3	AB 00
	1.29	AA 10 00 03	Set the output of the AlgoHub to sensor and algorithm.	AB 00

	1.30	AA 44 07 01 01 (CMD_DELAY = 500ms)	Enable Wearable Algorithm Suite AlgoHub	AB 00
READING SAMPLES REPORT IN OUTPUT FIFO	Host repeats the following commands every 1s			
	2.1	AA 14 00 [PPG1 PPG2 PPG3 PPG4 PPG5 PPG6 ACCLX ACCLY ACCLZ]  <b>25 [ppg, accel input frames] (max 25)</b> (CMD_DELAY = 5 ms)	Write LEDPPG, accel data to the algorithm input FIFO. PPG1: 3 bytes [LSB .. MSB] PPG2: 3 bytes [LSB .. MSB] PPG3: 3 bytes [LSB .. MSB] PPG4: 3 bytes [LSB .. MSB] PPG5: 3 bytes [LSB .. MSB] PPG6: 3 bytes [LSB .. MSB] Accelerometer set to ±8g: ACCLX: 2 byte accel X value (2's complement .001g): [LSB MSB] ACCLY: 2 byte accel Y value (2's complement .001g): [LSB MSB] ACCLZ: 2 byte accel Z value (2's complement .001g): [LSB MSB]	AB 00 mm nn  mmnn is the ppg, accel number of bytes received (16-bit unsigned)
		Delay for AlgoHub Processing – See “AlgoHub Response Duration”		
	2.2	AA 00 00	Read the AlgoHub status byte: Bit 0: Sensor comm error Bits 1 and 2: Reserved Bit 3: FIFO filled to threshold (DataRdyInt) Bit 4: Output FIFO overflow (FifoOutOvrInt) Bit 5: Input FIFO overflow (FifoInOverInt) Bit 6: AlgoHub busy (DevBusy) Bit 7: Reserved If DataRdyInt is set, proceed to the next step.	AB 00 08
	2.3	AA 12 00	Get the number of samples (nn) in the FIFO.	AB 00 nn
	2.4	AA 12 01 (CMD_DELAY = 5 ms)	Read the FIFO data for nn samples report. The format of the samples report is shown in the normal algorithm report table.	AB 00 nn*[Normal Algorithm Report]
		2.5, 2.6 are commanded if the AFE settings request flag is set		
	2.5	AA 47 07 27 (CMD_DELAY = 5 ms)	Read algorithm AFE settings requests.	[byte1 byte2 byte3 byte4 byte5]  byte1, byte2: bit 15: <b>0</b> : No LED current update request <b>1</b> : LED current update requested bits 14-0: LED current (15-bit unsigned, 0.1mA)  byte3: bit 7: <b>0</b> : No integration time update request <b>1</b> : Integration time update requested bits 6-0:

				AFE integration time: <b>0x00:</b> 14.8µs <b>0x01:</b> 29.4µs <b>0x02:</b> 58.7µs <b>0x03:</b> 117.3µs  byte4: bit 7: <b>0:</b> No sampling rate update request <b>1:</b> Sampling rate update requested bits 6-0: Sampling rate and averaging. <b>0x00:</b> 25sps, avg = 1 <b>0x01:</b> 50sps, avg = 2 <b>0x02:</b> 100sps, avg = 4 <b>0x03:</b> 200sps, avg = 8 <b>0x04:</b> 400sps, avg = 16  byte5: DAC offset bit 7: <b>0:</b> No DAC offset update request <b>1:</b> DAC offset update requested <b>0x00:</b> 0 uA <b>0x01:</b> 8 uA <b>0x02:</b> 16 uA <b>0x03:</b> 24 uA
	2.6	AA 47 07 28 (CMD_DELAY = 5 ms)	Clear algorithm AFE settings request flag. This command is sent after the host has honored the algorithm AFE requests.	AB 00 08
STOP	Host disables sensors and algorithm:			
	3.1	AA 44 07 00 01 (CMD_DELAY = 200ms)	Disable the WAS algorithm.	AB 00
	3.2	AA 46 26 (CMD_DELAY = 25ms)	Reset the AFE settings of the WAS algorithm.	AB 00

### AlgoHub AGC Normal Algorithm Report Configuration

The Automatic Gain Control (AGC) may be enabled by enabling AEC and disabling automatic PD calculation. The sequence of commands for enabling AGC is shown in the table below.

**Table 6. Host Commands—AlgoHub AGC and Read Normal Algorithm Report**



	#	HOST COMMAND (HEX)	COMMAND DESCRIPTION	RESPONSE (HEX)
START ALGORITHM		Host initializes the AlgoHub and starts the algorithm using following commands:		
	1.1	AA 02 00 (optional)	Read the operating mode	AB 00 00 application mode
	1.2	AA FF 03 (optional)	Read the AlgoHubSensorHub version for MAX32674C (x.y.z)	AB 00 32 YY ZZ
	1.3	AA 54 00	MAX32674C releases sensor SPI bus to Host. AlgoHub mode.	AB 00
	1.4	AA 10 01 01	Set the AlgoHub interrupt threshold.	AB 00
	1.5	AA 46 07 0B 01	Enable AEC	AB 00
	1.6	AA 46 07 12 00	Disable automatic PD calculation	AB 00
	1.7	AA 46 07 1A 00 03	Set PPG measurement 1 AFE integration time to 117.3us	AB 00
	1.8	AA 46 07 1A 01 03	Set PPG measurement 2 AFE integration time to 117.3us	AB 00
	1.9	AA 46 07 1A 02 03	Set PPG measurement 3 AFE integration time to 117.3us	AB 00
	1.10	AA 46 07 1B 00 01	Set PPG measurement 1 AFE sampling rate to 50sps, averaging to 2	AB 00
	1.11	AA 46 07 1B 01 04	Set PPG measurement 2 AFE sampling rate to 400sps, averaging to 16	AB 00
	1.12	AA 46 07 1B 02 04	Set PPG measurement 3 AFE sampling rate to 400sps, averaging to 16	AB 00
	1.13	AA 46 07 0F 00 7D	Set minimum PD current to 12.5uA	AB 00
	1.14	AA 46 07 10 01 38	Set initial PD current to 31.2uA	AB 00
	1.15	AA 46 07 0D 07 08	Set the adjusted target PD current period to 1800s	AB 00
	1.16	AA 46 07 0E 00 32	Set the HR motion magnitude threshold to 0.05g.	AB 00
	1.17	AA 46 07 23 00 00	Set the DAC offset for PD channel 1 for PPG measurement 1	AB 00
	1.18	AA 46 07 24 00 00	Set the DAC offset for PD channel 2 for PPG measurement 1	AB 00
	1.19	AA 46 07 1C 00 00	Set the maximum DAC offset for PPG measurement 1 to 0	AB 00
	1.20	AA 46 07 23 01 00	Set the DAC offset for PD channel 1 for PPG measurement 2	AB 00
	1.21	AA 46 07 24 01 00	Set the DAC offset for PD channel 2 for PPG measurement 2	AB 00
	1.22	AA 46 07 23 02 00	Set the DAC offset for PD channel 3 for PPG measurement 1	AB 00
	1.23	AA 46 07 24 02 00	Set the DAC offset for PD channel 2 for PPG measurement 3	AB 00
	1.24	AA 46 07 11 01 38	Set target PD current to 31.2uA	AB 00
	1.25	AA 46 07 0C 01	Enable SCD	AB 00
	1.26	AA 46 07 25 00 00 64	Set initial LED current to 10mA for PPG measurement 1	AB 00
	1.27	AA 46 07 25 01 00 C8	Set initial LED current to 20mA for PPG measurement 2	AB 00
	1.28	AA 46 07 25 02 00 C8	Set initial LED current to 20mA for PPG measurement 3	AB 00
	1.29	AA 10 00 03	Set the output of the AlgoHub to sensor and algorithm.	AB 00



	1.30	AA 44 07 01 01 (CMD_DELAY = 500ms)	Enable Wearable Algorithm Suite AlgoHub	AB 00
READING SAMPLES REPORT IN OUTPUT FIFO	Host repeats the following commands every 40ms			
	2.1	AA 14 00 [PPG1 PPG2 PPG3 PPG4 PPG5 PPG6 ACCLX ACCLY ACCLZ]  (CMD_DELAY = 16 ms)	Write LEDPPG, accel data to the algorithm input FIFO. PPG1: 3 bytes [MSB .. LSB] WHRM HR Channel 1 (green1) PPG2: 3 bytes [MSB .. LSB] WHRM HR Channel 2 (green2) PPG3: 3 bytes [MSB .. LSB] SpO2 IR Channel (IR) PPG4: 3 bytes [MSB .. LSB] SpO2 IR Channel (red) PPG5: 3 bytes [MSB .. LSB] PPG6: 3 bytes [MSB .. LSB] Accelerometer set to ±8g: ACCLX: 2 byte accel X value (2's complement .001g): [MSB LSB] ACCLY: 2 byte accel Y value (2's complement .001g): [MSB LSB] ACCLZ: 2 byte accel Z value (2's complement .001g): [MSB LSB]	AB 00 mm nn  mmnn is the ppg, accel number of bytes received (16-bit unsigned)
		Delay for 20 msec		
	2.2	AA 00 00	Read the AlgoHub status byte: Bit 0: Sensor comm error Bits 1 and 2: Reserved Bit 3: FIFO filled to threshold (DataRdyInt) Bit 4: Output FIFO overflow (FifoOutOvrInt) Bit 5: Input FIFO overflow (FifoInOverInt) Bit 6: AlgoHub busy (DevBusy) Bit 7: Reserved If DataRdyInt is set, proceed to the next step.	AB 00 08
	2.3	AA 12 01 (CMD_DELAY = 5 ms)	Read the data stored in the FIFO. The format of the samples report is shown in the normal algorithm report table.	AB 00 [Normal Algorithm Report]
	2.4, 2.5 are commanded if the AFE settings request flag is set			
	2.4	AA 47 07 27 (CMD_DELAY = 5 ms)	Read algorithm AFE settings requests.	[byte1 byte2 byte3 byte4 byte5]  byte1, byte2: bit 15: 0: No LED current update request 1: LED current update requested bits 14-0: LED current (15-bit unsigned, 0.1mA)  byte3: bit 7: 0: No integration time update request 1: Integration time update requested

				bits 6-0: AFE integration time: <b>0x00</b> : 14.8µs <b>0x01</b> : 29.4µs <b>0x02</b> : 58.7µs <b>0x03</b> : 117.3µs  byte4: bit 7: <b>0</b> : No sampling rate update request <b>1</b> : Sampling rate update requested bits 6-0: Sampling rate and averaging. <b>0x00</b> : 25sps, avg = 1 <b>0x01</b> : 50sps, avg = 2 <b>0x02</b> : 100sps, avg = 4 <b>0x03</b> : 200sps, avg = 8 <b>0x04</b> : 400sps, avg = 16  byte5: DAC offset bit 7: <b>0</b> : No DAC offset update request <b>1</b> : DAC offset update requested <b>0x00</b> : 0 uA <b>0x01</b> : 8 uA <b>0x02</b> : 16 uA <b>0x03</b> : 24 uA
	2.5	AA 47 07 28 (CMD_DELAY = 5 ms)	Clear algorithm AFE settings request flag. This command is sent after the host has honored the algorithm AFE requests.	AB 00 08
STOP	Host disables sensors and algorithm:			
	3.1	AA 44 07 00 01 (CMD_DELAY = 200ms)	Disable the WAS algorithm.	AB 00
	3.2	AA 46 26 (CMD_DELAY = 25ms)	Reset the AFE settings of the WAS algorithm.	AB 00

Configuration for AlgoHub Normal Algorithm Report when Host

### Configuration for AlgoHub Normal Algorithm Report when Host Implements AEC/AGC

If the system requirements does not need the AlgoHub to handle AEC/AGC and if the host firmware has implemented the AEC/AGC code, then the AlgoHub may disable the AlgoHub AEC (which also disables AlgoHub AGC) and configure the AlgoHub per the following table.

**Table 7. Host Commands—AlgoHub Normal Algorithm Report when Host Implements AEC/AGC**

	#	HOST COMMAND (HEX)	COMMAND DESCRIPTION	RESPONSE (HEX)
START ALGORITHM		Host initializes the AlgoHub and starts the algorithm using following commands:		
	1.1	AA 02 00 (optional)	Read the operating mode	AB 00 00 application mode
	1.2	AA FF 03 (optional)	Read the AlgoHubSensorHub version for MAX32674C (x.y.z)	AB 00 32 YY ZZ
	1.3	AA 54 00	MAX32674C releases sensor SPI bus to Host. AlgoHub mode.	AB 00
	1.4	AA 10 01 01	Set the AlgoHub interrupt threshold.	AB 00
	1.5	AA 46 07 0B 00	Enable AEC	AB 00
	1.6	AA 46 07 0C 01	Enable SCD	AB 00
	1.7	AA 46 07 25 00 00 64	Set initial LED current to 10mA for PPG measurement 1	AB 00
	1.8	AA 46 07 25 01 00 C8	Set initial LED current to 20mA for PPG measurement 2	AB 00
	1.9	AA 46 07 25 02 00 C8	Set initial LED current to 20mA for PPG measurement 3	AB 00
	1.10	AA 10 00 03	Set the output of the AlgoHub to sensor and algorithm.	AB 00
	1.11	AA 44 07 01 01 (CMD_DELAY = 500ms)	Enable Wearable Algorithm Suite AlgoHub	AB 00
READING SAMPLES REPORT IN OUTPUT FIFO		Host repeats the following commands every 40ms		
	2.1	AA 14 00 [PPG1 PPG2 PPG3 PPG4 PPG5 PPG6 ACCLX ACCLY ACCLZ]  (CMD_DELAY = 16 ms)	Write LEDPPG, accel data to the algorithm input FIFO. PPG1: 3 bytes [MSB .. LSB] WHRM HR Channel 1 (green1) (or if green is not used, then use (IR)) PPG2: 3 bytes [MSB .. LSB] WHRM HR Channel 2 (green2) (or if green is not used, then use (IR)) PPG3: 3 bytes [MSB .. LSB] SpO2 IR Channel (IR) PPG4: 3 bytes [MSB .. LSB] SpO2 IR Channel (red) PPG5: 3 bytes [MSB .. LSB] PPG6: 3 bytes [MSB .. LSB] Accelerometer set to ±8g: ACCLX: 2 byte accel X value (2's complement .001g): [MSB LSB] ACCLY: 2 byte accel Y value (2's complement .001g): [MSB LSB] ACCLZ: 2 byte accel Z value (2's complement .001g): [MSB LSB]	AB 00 mm nn  mmnn is the ppg, accel number of bytes received (16-bit unsigned)
		Delay for 20 msec		
	2.2	AA 00 00	Read the AlgoHub status byte: Bit 0: Sensor comm error Bits 1 and 2: Reserved Bit 3: FIFO filled to threshold (DataRdyInt) Bit 4: Output FIFO overflow (FifoOutOvrInt) Bit 5: Input FIFO overflow (FifoInOverInt) Bit 6: AlgoHub busy (DevBusy) Bit 7: Reserved	AB 00 08

			If DataRdyInt is set, proceed to the next step.	
2.3	AA 12 01 (CMD_DELAY = 5 ms)	Read the data stored in the FIFO. The format of the samples report is shown in the normal algorithm report table.	AB 00 [Normal Algorithm Report]	
2.4	AA 47 07 27 (CMD_DELAY = 5 ms)	Read algorithm AFE settings requests.	[byte1 byte2 byte3 byte4 byte5]  byte1, byte2: bit 15: <b>0</b> : No LED current update request <b>1</b> : LED current update requested bits 14-0: LED current (15-bit unsigned, 0.1mA)  byte3: bit 7: <b>0</b> : No integration time update request <b>1</b> : Integration time update requested bits 6-0: AFE integration time: <b>0x00</b> : 14.8µs <b>0x01</b> : 29.4µs <b>0x02</b> : 58.7µs <b>0x03</b> : 117.3µs  byte4: bit 7: <b>0</b> : No sampling rate update request <b>1</b> : Sampling rate update requested bits 6-0: Sampling rate and averaging. <b>0x00</b> : 25sps, avg = 1 <b>0x01</b> : 50sps, avg = 2 <b>0x02</b> : 100sps, avg = 4 <b>0x03</b> : 200sps, avg = 8 <b>0x04</b> : 400sps, avg = 16  byte5: DAC offset bit 7: <b>0</b> : No DAC offset update request <b>1</b> : DAC offset update requested <b>0x00</b> : 0 uA <b>0x01</b> : 8 uA	

				<b>0x02:</b> 16 uA <b>0x03:</b> 24 uA
	2.5	AA 47 07 28 (CMD_DELAY = 5 ms)	Clear algorithm AFE settings request flag. This command is sent after the host has honored the algorithm AFE requests.	AB 00 08
<b>STOP</b>	Host disables sensors and algorithm:			
	3.1	AA 44 07 00 01 (CMD_DELAY = 200ms)	Disable the WAS algorithm.	AB 00
	3.2	AA 46 26 (CMD_DELAY = 25ms)	Reset the AFE settings of the WAS algorithm.	AB 00

### AlgoHub Processing Response Duration

AlgoHub response duration depends on number of sample frames pushed to AlgoHub input FIFO in a single transaction. The CMD\_DELAY between the data feed command and the status byte is 5ms. The maximum number of sample frames that may be sent in one command is 25 per period. The AlgoHub processing time response time can be calculated as:

$$\text{Wait duration for AlgoHub Results} = 4\text{ms} + 2\text{ms} * (\# \text{ input frames in transaction})$$

For example, for 25 input samples, the wait time for 25 AlgoHub results is  $4 + 25*2 = 54\text{ms}$ .

## Configuring the SensorHub and Reading Processed Data

### I2C SensorHub Output FIFO Format

The SensorHub output FIFO may contain Sample Counter byte and/or Sensor Data, and/or WAS Algorithm Data and/or BPT Algorithm Data. The inclusion of these categories of output data is controlled by the “Set the output format of the sensor hub”, 0x10 0x00 command and the formats are listed below (The initial response byte, “Read Status Byte” is not shown).

**Table 8. SensorHub Sample Counter Byte Output FIFO Format**

DATA ITEM	# OF BYTES (MSB FIRST)	DESCRIPTION
Sample Counter byte	1	Counter which cycles from 0-0xFF.

**Table 9. SensorHub Sensor Data Output FIFO Format**

DATA SOURCE (24 bytes total)	DATA ITEM	# OF BYTES (MSB FIRST)	DESCRIPTION
Accelerometer (6 Bytes)	accelX	2	Two's complement. lsb = 0.001g
	accelY	2	Two's complement. lsb = 0.001g
	accelZ	2	Two's complement. lsb = 0.001g
MAX86176 PPG Data	PPG1 (PD1)	3	Green1

(18 Bytes)	PPG2 (PD1)	3	IR LED counts
	PPG3 (PD1)	3	Red LED counts
	PPG4 (PD2)	3	Green2
	PPG5 (PD2)	3	N/A
	PPG6 (PD2)	3	N/A

**Table 10. SensorHub WAS Algorithm Data Output FIFO Format**

DATA ITEM (20 bytes total)	# OF BYTE S (MSB FIRST)	DESCRIPTION
WAS algorithm operation mode (HR, SpO <sub>2</sub> )	1	WAS algorithm operation mode: 0: Continuous HRM and Continuous SpO <sub>2</sub> 1: Continuous HRM and One-Shot SpO <sub>2</sub> 2: Continuous HRM 3: Sampled HRM 4: Sampled HRM and One-Shot SpO <sub>2</sub> 5: Activity tracking 6: SpO <sub>2</sub> Calibration Data Collection
HR	2	10x Calculated heart rate
HR confidence	1	Calculated confidence level in %
RR	2	10x RR – inter-beat interval in ms Only shows a nonzero value when a new value is calculated.
RR confidence	1	Calculated confidence level of RR in % Only shows a nonzero value when a new value is calculated.
Activity class	1	Activity class: 0: Rest 1: Other 2: Walk 3: Run 4: Bike
R	2	1000x Calculated SpO <sub>2</sub> R value
SpO <sub>2</sub> confidence	1	Calculated SpO <sub>2</sub> confidence level in %
SpO <sub>2</sub>	2	10x Calculated SpO <sub>2</sub> %
SpO <sub>2</sub> percent complete	1	Calculation progress in percentage in one-shot mode of algorithm. In continuous mode, it is reported as zero and only jumps to 100 when the SpO <sub>2</sub> value is updated.
SpO <sub>2</sub> low signal quality flag	1	Shows the low quality of the PPG signal: 0: Good quality 1: Low quality
SpO <sub>2</sub> motion flag	1	Shows excessive motion: 0: No motion 1: Excessive motion
SpO <sub>2</sub> low PI flag	1	Shows the low perfusion index (PI) of the PPG signal: 0: Normal PI 1: Low PI
SpO <sub>2</sub> unreliable R flag	1	Shows the reliability of R: 0: Reliable 1: Unreliable
SpO <sub>2</sub> state	1	Reported status of the SpO <sub>2</sub> algorithm: 0: LED adjustment 1: Computation

		2: Success 3: Timeout
SCD state	1	Skin contact state: 0: Undetected 1: Off skin 2: On some subject 3: On skin

**Table 11. SensorHub WAS Extended Algorithm Data Output FIFO Format**

DATA ITEM (58 bytes total)	# OF BYTE S (MSB FIRST)	DESCRIPTION
Same as the previous table, WAS algorithm data	20	Same as the previous table, WAS algorithm data
Total walk steps	4	Total number of walking steps since the last reset
Total run steps	4	Total number of running steps since the last reset. (Applicable to wrist form factor only, MAX86141/0, MAXM86146)
Total energy exp in kcal	4	10x total energy expenditure since the last reset in kcal. (Applicable to wrist form factor only, MAX86141/0, MAXM86146)
Total AMR in kcal	4	10x total active energy expenditure since the last reset in kcal. (Applicable to wrist form factor only, MAX86141/0, MAXM86146)
AEC current update request for green1 algo channel	2	Bit15: is current update requested Bit0-14: requested current
AEC TINT request for green1 algo channel	1	Bit7: is TINT update requested Bit0-6: requested TINT
AEC sample average request for green1 algo channel	1	Bit7: is sampling average update requested Bit0-6: requested sampling average
AEC DAC offset request for green1 algo channel	1	Bit7: is DAC offset update requested Bit0-6: requested DAC offset
AEC current update request for green2 algo channel	2	Bit15: is current update requested Bit0-14: requested current
AEC TINT request for green2 algo channel	1	Bit7: is TINT update requested Bit0-6: requested TINT
AEC sample average request for green2 algo channel	1	Bit7: is sampling average update requested Bit0-6: requested sampling average
AEC DAC offset request for green2 algo channel	1	Bit7: is DAC offset update requested Bit0-6: requested DAC offset
AEC current update request for IR algo channel	2	Bit15: is current update requested Bit0-14: requested current



AEC TINT request for IR algo channel	1	Bit7: is TINT update requested Bit0-6: requested TINT
AEC sample average request for IR algo channel	1	Bit7: is sampling average update requested Bit0-6: requested sampling average
AEC DAC offset request for IR algo channel	1	Bit7: is DAC offset update requested Bit0-6: requested DAC offset
AEC current update request for RED algo channel	2	Bit15: is current update requested Bit0-14: requested current
AEC TINT request for RED algo channel	1	Bit7: is TINT update requested Bit0-6: requested TINT
AEC sample average request for RED algo channel	1	Bit7: is sampling average update requested Bit0-6: requested sampling average
AEC DAC offset request for RED algo channel	1	Bit7: is DAC offset update requested Bit0-6: requested DAC offset
AFE controller state for HR channels	1	Bit0-7: state
AFE controller output for high motion HR	1	Bit0: is high motion for high motion HR

**Table 12. SensorHub BPT Algorithm Data Output FIFO Format**

DATA ITEM (8 bytes total)	NUMBER OF BYTES (MSB FIRST)	DESCRIPTION
BPT status	1	BPT Status: 0: No signal 1: User calibration/estimation in progress 2: Success 3: Weak signal 4: Motion 5: Estimation failure 6: Calibration partially complete 7: Subject initialization failure 8: Initialization completed 9: Calibration reference BP trending error 10: Calibration reference Inconsistency 1 error 11: Calibration reference Inconsistency 2 error 12: Calibration reference Inconsistency 3 error 13: Calibration reference count mismatch 14: Calibration reference are out of limits (systolic 80 to 180, diastolic 50 to 120) 15: Number of calibrations exceed maximum 16: Pulse pressure out of range 17: Heart rate out of range

		18: Heart rate is above resting 19: Perfusion Index is out of range 20: Estimation error, try again 21: BPT estimate is out of range from calibration references (systolic +-30, diastolic +-20) 22: BPT estimate is beyond the maximum limits (systolic 80 to 180, diastolic 50 to 120)
BPT Progress	1	% complete
Systolic blood pressure	1	Estimated systolic blood pressure
Diastolic blood pressure	1	Estimated diastolic blood pressure
Bit-Mapped BPT Status	2	Bit-Mapped BPT Status 0x01: Pulse pressure out of range 0x02: Heart rate is above resting 0x04: Heart rate out of range 0x08: Perfusion Index is out of range 0x10: Estimation failure 0x20: : BPT estimate is out of range from calibration references (systolic +-30, diastolic +-20) 0x40: BPT estimate is beyond the maximum limits (systolic 80 to 180, diastolic 50 to 120)
Reserved	2	Reserved for future use

### SensorHub, Heart Rate, SpO2, AEC/AGC Report Configuration

Automatic Exposure Control (AEC) is Maxim's gain control algorithm that is superior to AGC. The AEC algorithm optimally maintains the best SNR range and power optimization. The targeted SNR range is maintained regardless of skin color or ambient temperature within the limits of the LED currents configurations; The AEC dynamically manages the appropriate register settings for sampling rate, LED current, pulse width and integration time.

In the example below, both the MAX32674C is configured to SensorHub mode with AEC, SCD enabled. The sequence of commands is shown in the table below.

**Table 13. Host Commands—SensorHub HR, SpO2, AEC/AGC Normal Algorithm Report**

	#	HOST COMMAND (HEX)	COMMAND DESCRIPTION	RESPONSE (HEX)
START ALGORITHM		Host initializes SensorHub and starts the algorithm using following commands:		
	0.1		Host sets SPI0_CSACC_N (P0.8) to high	
	0.2		Host initiates "Reset to Application Mode" sequence	
		Host initializes the SensorHub and starts the algorithm using following commands:		
	1.1	AA 02 00 (optional)	Read the operating mode	AB 00 00 application mode
	1.2	AA FF 03 (optional)	Read the AlgoHubSensorHub version for MAX32674C (x.y.z)	AB 00 32 YY ZZ
	1.3	AA 10 01 01	Set the SensorHub interrupt threshold.	AB 00
	1.4	AA 54 01	MAX32674C releases sensor SPI bus to MAX32674C. SensorHub mode.	AB 00
	1.6	AA 10 02 25	Set the samples report period to 1s	AB 00
		AA 10 02 1	Set the samples report period to 40ms	
	1.6	AA 40 06 1E 1F (optional)	Write 0x1F to MAX86176 register 0x1E	AB 00
	1.7	AA 40 06 1D 05 (optional)	Write 0x05 to MAX86176 register 0x1D	AB 00
	1.8	AA 40 06 11 07 (optional)	Write 0x07 to MAX86176 register 0x11	AB 00
	1.9	AA 40 06 12 04 (optional)	Write 0x04 to MAX86176 register 0x12	AB 00
	1.10	AA 40 06 90 02 (optional)	Write 0x90 to MAX86176 register 0x90	AB 00
	1.11	AA 50 08 0B 01	Enable AEC	AB 00
	1.12	AA 50 08 12 00 (alternate choice)	Disable automatic PD calculation for AGC	AB 00
		AA 50 08 12 01	Enable automatic PD calculation for AEC	
	1.13	AA 50 08 1A 00 03 (optional)	Set PPG measurement 1 AFE integration time to 117.3us	AB 00
	1.14	AA 50 08 1A 01 03 (optional)	Set PPG measurement 2 AFE integration time to 117.3us	AB 00
	1.15	AA 50 08 1A 02 03 (optional)	Set PPG measurement 3 AFE integration time to 117.3us	AB 00
	1.16	AA 50 08 1B 00 01 (optional)	Set PPG measurement 1 AFE sampling rate to 50sps, averaging to 2	AB 00
	1.17	AA 50 08 1B 01 04 (optional)	Set PPG measurement 2 AFE sampling rate to 400sps, averaging to 16	AB 00
	1.18	AA 50 08 1B 02 04 (optional)	Set PPG measurement 3 AFE sampling rate to 400sps, averaging to 16	AB 00
	1.19	AA 50 08 0F 00 7D (optional)	Set minimum PD current to 12.5uA	AB 00
	1.20	AA 50 08 10 01 38 (optional)	Set initial PD current to 31.2uA	AB 00
	1.21	AA 50 08 0D 07 08 (optional)	Set the adjusted target PD current period to 1800s	AB 00
	1.22	AA 50 08 0E 00 32 (optional)	Set the HR motion magnitude threshold to 0.05g.	AB 00

1.23	AA 50 08 23 00 00 (optional)	Set the DAC offset for PD channel 1 for PPG measurement 1	AB 00
1.24	AA 50 08 24 00 00 (optional)	Set the DAC offset for PD channel 2 for PPG measurement 1	AB 00
1.25	AA 50 08 1C 00 00 (optional)	Set the maximum DAC offset for PPG measurement 1 to 0	AB 00
1.26	AA 50 08 23 01 00 (optional)	Set the DAC offset for PD channel 1 for PPG measurement 2	AB 00
1.27	AA 50 08 24 01 00 (optional)	Set the DAC offset for PD channel 2 for PPG measurement 2	AB 00
1.28	AA 50 08 23 02 00 (optional)	Set the DAC offset for PD channel 3 for PPG measurement 1	AB 00
1.29	AA 50 08 24 02 00 (optional)	Set the DAC offset for PD channel 2 for PPG measurement 3	AB 00
1.30	AA 50 08 25 00 00 64 (optional)	Set initial LED current to 10mA for PPG measurement 1	AB 00
1.31	AA 50 08 25 01 00 C8 (optional)	Set initial LED current to 20mA for PPG measurement 2	AB 00
1.32	AA 50 08 25 02 00 C8 (optional)	Set initial LED current to 20mA for PPG measurement 3	AB 00
1.33	AA 50 08 11 01 38 (optional)	Set target PD current to 31.2uA	AB 00
1.34	AA 50 08 0C 01 (optional default)	Enable SCD	AB 00
1.35	AA 10 00 03 (alternate choice)	Set the output format to Sensor Data and Algorithm	AB 00
	AA 10 00 07	Set the output format to Sample Counter byte, Sensor Data and Algorithm	
1.36	AA 40 06 10 01 (optional)	Write 0x01 to MAX86176 register 0x10	AB 00
1.37	AA 40 06 10 02 (optional)	Write 0x02 to MAX86176 register 0x10	AB 00
1.38	AA 40 06 15 1E (optional)	Write 0x15 to MAX86176 register 0x1E	AB 00
1.39	AA 40 06 18 80 (optional)	Write 0x18 to MAX86176 register 0x80	AB 00
1.40	AA 44 04 01 00 (CMD_DELAY = 50ms)	Enable the accelerometer.	AB 00
1.41	AA 44 06 01 00 (CMD_DELAY = 500ms)	Enable MAX86176 AFE.	AB 00
1.42	AA 50 08 40 01 (optional default)	Set the biometric operation mode to WAS	AB 00
1.43	AA 50 08 0A 00 (optional default)	Set the WAS algorithm operation mode to Continuous HRM + Continuous SpO2	AB 00
1.44	AA 52 08 02 (CMD_DELAY = 500ms) (alternate choice)	Write enable/disable the Wearable Algorithm Suite (WHRM+WSpO2) extended algorithm	AB 00
	AA 52 08 01 (CMD_DELAY = 500ms)	Write enable/disable the Wearable Algorithm Suite (WHRM+WSpO2) algorithm	

	1.45	AA 41 00 FF (optional)	Read register FF (PART_ID) of MAX86176	AB 00 [39]
	1.46	AA 41 04 0F (accel data is required LIS2DS12I)	Read register 0F (WHO_AM_I) of LIS2DS12	AB 00 [43]
READING SAMPLES REPORT IN OUTPUT	Host repeats the following commands every 40ms			
	2.1	AA 00 00	Read the sensor hub status byte: Bit 0: Sensor comm error Bits 1 and 2: Reserved Bit 3: FIFO filled to threshold (DataRdyInt) Bit 4: Output FIFO overflow (FifoOutOvrInt) Bit 5: Input FIFO overflow (FifoInOverInt) If DataRdyInt is set, proceed to the next step.	AB 00 08
	2.2	AA 12 01	Read the data stored in the FIFO	AB 00 [ sample counter byte] [sensor data] [algorithm data]
STOP	Host disables sensors and algorithm:			
	3.1	AA 44 04 00 (CMD_DELAY = 50ms)	Disable the accelerometer.	AB 00
	3.2	AA 44 06 00 (CMD_DELAY = 200ms)	Disable the MAX86176 AFE.	AB 00
	3.2	AA 52 08 00 (CMD_DELAY = 200ms)	Disable the Wearable Algorithm Suite (WHRM+WSpO2) algorithm	AB 00

### SensorHub, PPG Raw Report Configuration

The host may configure the SensorHub to output raw data (no algorithm) by enabling the AFE. The table below lists the set of commands that are needed to obtain the PPG raw data.

**Table 14. Host Commands— SensorHub PPG Raw Data Report**

	#	HOST COMMAND (HEX)	COMMAND DESCRIPTION	RESPONSE (HEX)
START ALGORITHM	Host initializes the SensorHub and AFE for raw data report:			
	1.1	AA 02 00 (optional)	Read the operating mode	AB 00 00 application mode
	1.2	AA FF 03 (optional)	Read the AlgoHubSensorHub version for MAX32674C (x.y.z)	AB 00 32 YY ZZ
	1.3	AA 54 01	MAX32674C releases sensor SPI bus to MAX32674C. SensorHub mode.	AB 00
	1.4	AA 10 00 03	Set the output format to Raw Data	AB 00
		AA 10 00 07	Set the output format to Sample Counter byte, Raw Data	
	1.5	AA 10 01 01	Set the output format to Sample Counter byte, Sensor Data and Algorithm	AB 00
	1.6	AA 10 02 01	set report period for Sensor Hub to 1	AB 00
	1.7	AA 44 04 01 00	Enable the accelerometer.	AB 00

		(CMD_DELAY = 50ms)		
	1.8	AA 44 06 01 00 (CMD_DELAY = 500ms)	Enable MAX86176 AFE.	AB 00
	1.9	AA 50 08 40 00	Set algo class operation mode to RAW mode	AB 00
	1.10	AA 52 08 01  (CMD_DELAY = 500ms)	Write enable the Wearable Algorithm Suite  Note: In RAW mode Algo class will work in bypass fashion.	AB 00
	1.11	AA 40 06 XX YY	Write desired setting to MEAS settings of channels for MEAS1 to MEAS3 for samp average, integration time, Led current, Dac offset.  <b>Note: do not change</b> 1. Firing sequence 2. Sampling Frequency	AB 00
READING SAMPLES REPORT IN OUTPUT FIFO	Host repeats the following commands every N ms			
	2.1	AA 00 00	Read the sensor hub status byte: Bit 0: Sensor comm error Bits 1 and 2: Reserved Bit 3: FIFO filled to threshold (DataRdyInt) Bit 4: Output FIFO overflow (FifoOutOvrInt) Bit 5: Input FIFO overflow (FifoInOverInt) If DataRdyInt is set, proceed to the next step.	AB 00 08
	2.2	AA 12 00	Read the number of samples available.	AB 00 01
	2.3	AA 12 01	Read the data stored in the FIFO	AB 00 ([sample counter byte] [raw data])x(number of samples)
STOP	Host disables sensors and algorithm:			
	3.1	AA 44 06 00 (CMD_DELAY = 200ms)	Disable the MAX86176 AFE.	AB 00

The SensorHub PPG raw data output FIFO format is listed below. The “Type” is defined in Table 12 of the MAX86176 Datasheet.

**Table 15. SensorHub PPG Raw Data Output FIFO Format**

DATA SOURCE (24 bytes)	DATA ITEM	# OF BYTES (MSB FIRST)	DESCRIPTION
MAX86176 PPG Data (18 Bytes)	Green (PD1)	3	[23:20]Type [19:0]: Green1
	Green (PD2)	3	[23:20]Type [19:0]: Green2
	IR (PD1)	3	[23:20]Type [19:0]: IR1
	IR (PD2)	3	[23:20]Type [19:0]: IR2
	Red (PD1)	3	[23:20]Type [19:0]: Red1
	Red (PD2)	3	[23:20]Type [19:0]: Red2

### SensorHub BPT, Heart Rate, SpO2 Report Configuration



### BPT User Calibration Vector Report

A multi-point (maximum of six) calibration process is required for each user to tune the BPT estimation algorithm. This step is done by first measuring user's reference systolic and diastolic BP using a medically approved device. Next, the time and date, the cal\_index, the cuff BP reference values are provided to the MAX32674C before each of the user calibration procedures are performed. This process takes about 1 minute. After the user calibration procedure is complete (BTP progress is 100% and BPT status is 2), a calibration vector is generated. For optimal results, it is recommended that five sets of these vectors (cal\_index 0 to 4) be generated for the user. The user calibration vector is valid if the MAX32674C has not re-started, and the user has not changed. However, after each calibration, it is recommended that the host (and app) read each calibration vector and store it in flash (or to a file) for a user and reload it to the MAX32674C prior to future estimation measurements.

The user calibration procedure is needed once a month to ensure the accuracy of BP estimation for the current user. The table below shows the sequence of commands for completing BPT user calibration.

**Table 16. Host Commands—BPT User Calibration Vector Report**

	#	HOST COMMAND (HEX)	COMMAND DESCRIPTION	RESPONSE (HEX)
INITIALIZE AND START BPT USER CALIBRATION	1.1	AA 02 00 (optional)	Read the operating mode	AB 00 00 application mode
	1.2	AA FF 03 (optional)	Read the AlgoHubSensorHub version for MAX32674C (x.y.z)	AB 00 32 YY ZZ
	1.3	AA 54 01	Set MAX32674C SPI bus to sensor hub mode.	AB 00
	1.4	AA 50 08 62 00 75 4C (example cal_index 0, reference systolic 117 and diastolic 76) (CMD_DELAY is 5 ms)	Set the cal_index, reference systolic, reference diastolic.	AB 00
	1.5	AA 50 08 0B 01	Enable AEC	AB 00
	1.6	AA 50 08 12 01	Enable automatic calculation of target PD current	AB 00
	1.7	AA 50 08 1A 00 03 (optional)	Set PPG measurement 1 initial integration time to 117.3us	AB 00
	1.8	AA 50 08 15 01 03 (optional)	Set PPG measurement 2 maximum integration time to 117.3us	AB 00
	1.9	AA 50 08 15 02 03 (optional)	Set PPG measurement 3 maximum integration time to 117.3us	AB 00
	1.10	AA 50 08 1B 00 01 (optional)	Set PPG measurement 1 initial sampling rate to 1	AB 00
	1.11	AA 50 08 1B 01 04 (optional)	Set PPG measurement 2 initial sampling rate to 4	AB 00
	1.12	AA 50 08 14 01 04 (optional)	Set PPG measurement 2 minimum sampling rate to 4	AB 00
	1.13	AA 50 08 16 01 10 (optional)	Set PPG measurement 2 maximum sampling rate to 16	AB 00
	1.14	AA 50 08 1B 02 04 (optional)	Set PPG measurement 3 initial sampling rate to 4	
	1.15	AA 50 08 14 02 04 (optional)	Set PPG measurement 3 minimum sampling rate to 4	AB 00



1.16	AA 50 08 16 02 10 (optional)	Set PPG measurement 3 maximum sampling rate to 16	
1.17	AA 50 08 0F 00 7D (recommended)	Set minimum PD current to 12.5uA	AB 00
1.18	AA 50 08 10 01 38 (recommended)	Set initial PD current to 31.2uA	AB 00
1.19	AA 50 08 0D 07 08 (optional)	Set the adjusted target PD current period to 1800s	AB 00
1.20	AA 50 08 0E 00 32 (optional)	Set the HR motion magnitude threshold to 0.05g.	AB 00
1.21	AA 50 08 23 00 00 (optional)	Set the DAC offset for PD channel 1 for PPG measurement 1	AB 00
1.22	AA 50 08 24 00 00 (optional)	Set the DAC offset for PD channel 2 for PPG measurement 1	AB 00
1.23	AA 50 08 1C 00 00 (optional)	Set the maximum DAC offset for PPG measurement 1 to 0	AB 00
1.24	AA 50 08 23 01 00 (optional)	Set the DAC offset for PD channel 1 for PPG measurement 2	AB 00
1.25	AA 50 08 24 01 00 (optional)	Set the DAC offset for PD channel 2 for PPG measurement 2	AB 00
1.26	AA 50 08 23 02 00 (optional)	Set the DAC offset for PD channel 3 for PPG measurement 1	AB 00
1.27	AA 50 08 24 02 00 (optional)	Set the DAC offset for PD channel 2 for PPG measurement 3	AB 00
1.28	AA 50 08 25 00 00 64 (optional)	Set initial LED current to 10mA for PPG measurement 1	AB 00
1.29	AA 50 08 25 01 00 C8 (optional)	Set initial LED current to 20mA for PPG measurement 2	AB 00
1.30	AA 50 08 25 02 00 C8 (optional)	Set initial LED current to 20mA for PPG measurement 3	AB 00
1.31	AA 50 08 11 01 38	Set target PD current to 31.2uA	AB 00
1.32	AA 50 08 0C 01	Enable SCD	AB 00
1.33	AA 10 00 07	Set the output format to Sample Counter byte, Sensor Data and Algorithm	AB 00
1.34	AA 10 02 05	Set the samples report period to 40ms (minimum is 32ms for BPT).	AB 00
1.35	AA 40 06 10 01 (optional)	Write 0x01 to MAX86176 register 0x10	AB 00
1.36	AA 40 06 10 02 (optional)	Write 0x02 to MAX86176 register 0x10	AB 00
1.37	AA 40 06 15 1E (optional)	Write 0x15 to MAX86176 register 0x1E	AB 00
1.38	AA 40 06 18 80 (optional)	Write 0x18 to MAX86176 register 0x80	AB 00
1.39	AA 40 06 19 9F (optional)	Write 0x19 to MAX86176 register 0x9F	AB 00
1.40	AA 40 06 1A 3F (optional)	Write 0x1A to MAX86176 register 0x3F	AB 00
1.41	AA 40 06 1C 20 (optional)	Write 0x1C to MAX86176 register 0x20	AB 00
1.42	AA 40 06 1D 05 (optional)	Write 0x1D to MAX86176 register 0x05	AB 00
1.43	AA 40 06 1E 1F (optional)	Write 0x1E to MAX86176 register 0x1F	AB 00
1.44	AA 40 06 0D E8 (optional)	Write 0x0D to MAX86176 register 0xE8	AB 00
1.45	AA 40 06 86 80 (optional)	Write 0x86 to MAX86176 register 0x80	AB 00
1.46	AA 44 04 01 00 (CMD_DELAY = 50ms)	Enable the accelerometer.	AB 00
1.47	AA 44 06 01 00 (CMD_DELAY = 500ms)	Enable MAX86176 AFE.	AB 00
1.48	AA 50 08 40 03	Set the biometric operation mode to WAS and BPT	AB 00

	1.49	AA 50 08 51 00	Set the BPT run mode to calibration	AB 00
	1.50	52 08 01 (CMD_DELAY = 500ms) (alternative)	Enable the biometric operation mode for BPT+WAS report algorithm	AB 00
		52 08 02 (CMD_DELAY = 500ms)	Enable the biometric operation mode for BPT+WAS extended report algorithm	
	1.51	AA 41 00 FF (optional)	Read register FF (PART_ID) of MAX86176	AB 00 [39]
	1.52	AA 41 04 0F (accel data is required LIS2DS12)	Read register 0F (WHO_AM_I) of LIS2DS12	AB 00 [43]
READING SAMPLES REPORT IN OUTPUT FIFO	Host repeats the following commands every 40ms			
	2.1	AA 00 00	Read the sensor hub status byte: Bit 0: Sensor comm error Bits 1 and 2: Reserved Bit 3: FIFO filled to threshold (DataRdyInt) Bit 4: Output FIFO overflow (FifoOutOvrInt) Bit 5: Input FIFO overflow (FifoInOverInt) If DataRdyInt is set, proceed to the next step.	AB 00 08 (Bit 3, DataRdyInt set)
	2.2	AA 12 01	Read the data stored in the FIFO	AB 00 [sample counter byte] [sensor data] [WAS algorithm extended data] [BPT algorithm data]
	Repeat 2.1 to 2.2 until BPT status is 2 and BPT progress is at 100%. If there's an error status or host timeout, then host should ask the user to restart the user calibration for this cal_index.			
STOP	Host disables sensors and algorithm:			
	3.1	AA 44 04 00 (CMD_DELAY = 50ms)	Disable the accelerometer.	AB 00
	3.2	AA 44 06 00 (CMD_DELAY = 200ms)	Disable the MAX86176 AFE.	AB 00
	3.3	AA 52 08 00 (CMD_DELAY = 200ms)	Disable the algorithm	AB 00
SAVE	4.1	AA 51 08 60	Read user calibration vector. Host saves the calibration vectors to user_name_profile[cal_index].	AB 00 [240 bytes of user calibration vector]
Repeat 1.1-1.50, 2.1-2.2, 3.1-3.3, 4.1 four more times (with cal_index increasing by one each time), if the user is performing five calibrations.				

#### **BPT Estimation Report and Heart Rate, SpO2 Report**

After the user calibration procedure is completed or restored, the algorithm can provide BPT, heart rate and SpO2 estimates using the configuration commands in the table below.

**Table 17. Host Commands—BPT Estimation Report and Heart Rate, SpO2 Report**

	#	HOST COMMAND (HEX)	COMMAND DESCRIPTION	RESPONSE (HEX)
17	1.1	AA 02 00 (optional)	Read the operating mode	AB 00 00

			application mode
1.2	AA FF 03 (optional)	Read the AlgoHubSensorHub version for MAX32674C (x.y.z)	AB 00 32 YY ZZ
1.3	AA 54 01	Set MAX32674C SPI bus to sensor hub mode.	AB 00
1.4	AA 50 08 0B 01	Enable AEC	AB 00
1.5	AA 50 08 12 00 (alternative)	Disable automatic calculation of target PD current for AGC	AB 00
	AA 50 08 12 01	Enable automatic calculation of target PD current for AEC	
1.6	AA 50 08 1A 00 03 (optional)	Set PPG measurement 1 initial integration time to 117.3us	AB 00
1.7	AA 50 08 15 01 03 (optional)	Set PPG measurement 2 maximum integration time to 117.3us	AB 00
1.8	AA 50 08 15 02 03 (optional)	Set PPG measurement 3 maximum integration time to 117.3us	AB 00
1.9	AA 50 08 1B 00 01 (optional)	Set PPG measurement 1 initial sampling rate to 1	AB 00
1.10	AA 50 08 1B 01 04 (optional)	Set PPG measurement 2 initial sampling rate to 4	AB 00
1.11	AA 50 08 14 01 04 (optional)	Set PPG measurement 2 minimum sampling rate to 4	AB 00
1.12	AA 50 08 16 01 10 (optional)	Set PPG measurement 2 maximum sampling rate to 16	AB 00
1.13	AA 50 08 1B 02 04 (optional)	Set PPG measurement 3 initial sampling rate to 4	
1.14	AA 50 08 14 02 04 (optional)	Set PPG measurement 3 minimum sampling rate to 4	AB 00
1.15	AA 50 08 16 02 10 (optional)	Set PPG measurement 3 maximum sampling rate to 16	
1.16	AA 50 08 0F 00 7D (recommended)	Set minimum PD current to 12.5uA	AB 00
1.17	AA 50 08 10 01 38 (recommended)	Set initial PD current to 31.2uA	AB 00
1.18	AA 50 08 0D 07 08 (optional)	Set the adjusted target PD current period to 1800s	AB 00
1.19	AA 50 08 0E 00 32 (optional)	Set the HR motion magnitude threshold to 0.05g.	AB 00
1.20	AA 50 08 23 00 00 (optional)	Set the DAC offset for PD channel 1 for PPG measurement 1	AB 00
1.21	AA 50 08 24 00 00 (optional)	Set the DAC offset for PD channel 2 for PPG measurement 1	AB 00
1.22	AA 50 08 1C 00 00 (optional)	Set the maximum DAC offset for PPG measurement 1 to 0	AB 00
1.23	AA 50 08 23 01 00 (optional)	Set the DAC offset for PD channel 1 for PPG measurement 2	AB 00
1.24	AA 50 08 24 01 00 (optional)	Set the DAC offset for PD channel 2 for PPG measurement 2	AB 00
1.25	AA 50 08 23 02 00 (optional)	Set the DAC offset for PD channel 3 for PPG measurement 1	AB 00
1.26	AA 50 08 24 02 00 (optional)	Set the DAC offset for PD channel 2 for PPG measurement 3	AB 00

	1.27	AA 50 08 25 00 00 64 (optional)	Set initial LED current to 10mA for PPG measurement 1	AB 00
	1.28	AA 50 08 25 01 00 C8 (optional)	Set initial LED current to 20mA for PPG measurement 2	AB 00
	1.29	AA 50 08 25 02 00 C8 (optional)	Set initial LED current to 20mA for PPG measurement 3	AB 00
	1.30	AA 50 08 11 01 38 (recommended)	Set target PD current to 31.2uA	AB 00
	1.31	AA 50 08 0C 01	Enable SCD	AB 00
	1.31	AA 10 00 03 (alternative)	Set the output format to Sensor Data and Algorithm	AB 00
	1.32		Set the output format to Sample Counter byte, Sensor Data and Algorithm	
	1.33	AA 10 00 07	Set the samples report period to 4 (minimum of 4 for BPT).	AB 00
		AA 10 02 04		
	1.34	AA 40 06 10 01 (optional)	Write 0x01 to MAX86176 register 0x10	AB 00
	1.35	AA 40 06 10 02 (optional)	Write 0x02 to MAX86176 register 0x10	AB 00
	1.36	AA 40 06 15 1E (optional)	Write 0x15 to MAX86176 register 0x1E	AB 00
	1.37	AA 40 06 18 80 (optional)	Write 0x18 to MAX86176 register 0x80	AB 00
	1.38	AA 40 06 19 9F (optional)	Write 0x19 to MAX86176 register 0x9F	AB 00
	1.39	AA 40 06 1A 3F (optional)	Write 0x1A to MAX86176 register 0x3F	AB 00
	1.40	AA 40 06 1C 20 (optional)	Write 0x1C to MAX86176 register 0x20	AB 00
	1.41	AA 40 06 1D 05 (optional)	Write 0x1D to MAX86176 register 0x05	AB 00
	1.42	AA 40 06 1E 1F (optional)	Write 0x1E to MAX86176 register 0x1F	AB 00
	1.43	AA 40 06 0D E8 (optional)	Write 0x0D to MAX86176 register 0xE8	AB 00
	1.44	AA 40 06 86 80 (optional)	Write 0x86 to MAX86176 register 0x80	AB 00
	1.45	AA 44 04 01 00 (CMD_DELAY = 50ms)	Enable the accelerometer.	AB 00
	1.46	AA 44 06 01 00 (CMD_DELAY = 500ms)	Enable MAX86176 AFE.	AB 00
	1.47	AA 50 08 40 03	Set the biometric operation mode to WAS and BPT	AB 00
	1.48	AA 50 08 51 01	Set the BPT run mode to estimation	AB 00
	1.49	52 08 01 (CMD_DELAY = 500ms) (alternative)	Enable the biometric operation mode for BPT+WAS report algorithm	AB 00
	52 08 02 (CMD_DELAY = 500ms)	Enable the biometric operation mode for BPT+WAS extended report algorithm		
READING SAMPLES REPORT IN OUTPUT FIFO	Host repeats the following commands every 40ms			
	2.1	AA 00 00	Read the sensor hub status byte: Bit 0: Sensor comm error Bits 1 and 2: Reserved Bit 3: FIFO filled to threshold (DataRdyInt) Bit 4: Output FIFO overflow (FifoOutOvrInt) Bit 5: Input FIFO overflow (FifoInOverInt) If DataRdyInt is set, proceed to the next step.	AB 00 08 (Bit 3, DataRdyInt set)
	2.2	AA 12 01	Read the data stored in the FIFO	AB 00 [sample counter byte] [sensor data] [WAS algorithm extended data]

				[BPT algorithm data]
		Repeat 2.1 to 2.2. BPT estimates are valid if BPT status is 2 and BPT progress is at 100%.		
STOP		Host disables sensors and algorithm:		
	3.1	AA 44 04 00 (CMD_DELAY = 50ms)	Disable the accelerometer.	AB 00
	3.2	AA 44 06 00 (CMD_DELAY = 200ms)	Disable the MAX86176 AFE.	AB 00
	3.3	AA 52 08 00 (CMD_DELAY = 200ms)	Disable the algorithm	AB 00
SAVE	4.1	AA 51 08 60	Read user calibration vector. Host saves the calibration vectors to user_name_profile[cal_index].	AB 00 [240 bytes of user calibration vector]

#### **BPT Restore User Calibration Vectors**

If the MAX32674C or app has been reset, or if the user\_name has changed in the app, then the BPT user calibration vectors need to be restored to the MAX32674C. The user may run the BPT estimation after these vectors have been restored. The table below illustrates how to restore the user calibration vectors for a given user that has five calibration vectors.

**Table 18. Host Commands—BPT Restore User Calibration Vectors**

	#	HOST COMMAND (HEX)	COMMAND DESCRIPTION	RESPONSE (HEX)
BPT RESTORE USER CALIBRATION VECTORS	1.1	AA 50 08 62 00	Set the the cal_index to 0.	AB 00
	1.2	AA 50 08 60 [240 bytes]	Set the BPT user calibration vector (using cal_index 0).	AB 00
	1.3	AA 50 08 62 01	Set the the cal_index to 1.	AB 00
	1.4	AA 50 08 60 [240 bytes]	Set the BPT user calibration vector (using cal_index 0).	AB 00
	1.4	AA 50 08 62 02	Set the the cal_index to 2.	AB 00
	1.5	AA 50 08 60 [240 bytes]	Set the BPT user calibration vector (using cal_index 0).	AB 00
	1.6	AA 50 08 62 03	Set the the cal_index to 3.	AB 00
	1.7	AA 50 08 60 [240 bytes]	Set the BPT user calibration vector (using cal_index 0).	AB 00
	1.8	AA 50 08 62 04	Set the the cal_index to 4.	AB 00

## Configuring a New MAX32674C

The MAX32674C is not pre-programmed with the application algorithm. See the section “Application .msbl Programming Sequence” to flash the latest .msbl application algorithm to the MAX32674.

The default settings for the MFIO pin must be configured once before the .msbl is flashed. The table below configures and saves the MFIO pin to P0.0.


**Table 19. Configure the MFIO pin and Save the Configuration**

	#	HOST COMMAND (HEX)	COMMAND DESCRIPTION	RESPONSE (HEX)
CONFIGURE THE MFIO PIN		Host must initialize the factory fresh MAX32674C once.		
	1.0		Reset to bootloader mode	
	1.1	AA 01 08	Stay in bootloader mode. This command must be sent within 20ms after reset.	AB AA
	1.2	AA 02 00 (optional)	Read the operating mode	AB AA 08 bootloader mode
	1.3	AA 82 01 00 01	Configure bootloader to use the MFIO pin to determine if the device enters bootloader or application mode after a reset.	AB AA
	1.4	AA 82 01 01 00	Configure the MFIO pin to P0.0	AB AA
	1.5	AA 82 01 02 00	Configure the MFIO pin to active low for bootloader mode.	AB AA
	1.6	AA 82 00 (CMD_DELAY = 300ms)	Save the configuration	AB AA
	1.7	Reset (see the next section on how to reset)	The MAX32674C must be reset to store the configure	

## Application .msbl Programming Sequence

To program the MAX32674C application .msbl, the host microprocessor places the MAX32674C into bootloader mode and the host may implement the software to flash the .msbl file (I2C, MFIO, RSTN lines are connected between the host and the MAX32674C). The MAX32674C uses the 8-bit slave address of 0xAA. Each page sent includes 16 CRC bytes for that page, so there are 8208 bytes per page sent in the payload of the message. The number of pages is located at address 0x44 in the .msbl file. Values for the number of pages, initialization vector, authorization bytes and page contents, might be different for the latest .msbl, but the locations of these values in the .msbl file remain the same; The example below used the data from the v50.1.0 .msbl file. There are additional bytes in the .msbl past the last page; these are the file checksum bytes. Since the bootloader uses the commands listed below and it does not accept files, the file checksum bytes are not used by the bootloader.

**Table 20. Annotated I2C Trace for Flashing the Application**

HOST COMMAND	COMMAND DESCRIPTION	READ MAX32674C RESPONSE	RESPONSE DESCRIPTION
Sequence the MAX32674C to enter bootloader mode. *			
 <p>Figure 1. Sequence to enter bootloader mode.</p>			



0xAA 0x02 0x00†	Read mode.	0xAB 0xAA 0x08	No error. Mode is bootloader.
0xAA 0x81 0x01†	Read bootloader page size.	0xAB 0xAA 0x20 0x00	No error. Page size is 8192.
0xAA 0x80 0x02 0x00 0x21*	Bootloader flash. Set the “number of pages” to 33 based on the value at byte 0x44 from the application .msbl file.	0xAB 0xAA	No error.

00000044 ba 46 34 2b 21 00 00 20 04 00 00 00 c7 f4 20 8b

Figure 2. Page number byte 0x44 from the .msbl file.

0xAA 0x80 0x00 0x8E 0xA2 0x9D 0x1A 0xE2 0x8F 0x7F 0x25 0x5E 0x0B 0x91*	Bootloader flash. Set the initialization vector bytes to the 0x28 to 0x32 values from the .msbl file.	0xAB 0xAA	No error.
---	---	-----------	-----------

00000020 00 00 00 00 00 00 00 00 8e a2 9d 1a e2 8f 7f 25  
00000032 5e 0b 91 00 0d e8 f8 12 7e 2e 8e d8 a9 a3 f1 60

Figure 3. Initialization vector bytes 0x28 to 0x32 from the .msbl file.

0xAA 0x80 0x01 0x0D 0xE8 0xF8 0x12 0x7E 0x2E 0x8E 0xD8 0xA9 0xA3 0xF1 0x60 0xBA 0x46 0x34 0x2B*	Bootloader flash. Set the authentication bytes to the 0x34 to 0x43 values from the .msbl file.	0xAB 0xAA	No error.
---	--	-----------	-----------

00000030 5e 0b 91 00 0d e8 f8 12 7e 2e 8e d8 a9 a3 f1 60  
00000043 ba 46 34 2b 21 00 00 20 04 00 00 00 c7 f4 20 8b

Figure 4. Authentication bytes 0x34 to 0x43 from the .msbl file.

0xAA 0x80 0x03* (CMD_DELAY = 1400ms)	Bootloader flash. Erase application.	0xAB 0xAA	No error.
0xAA 0x80 0x04 0xC7 0xF4 0x20 ... 0xD4 0x75 0xBF* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0x4C to 0x205B from the .msbl file (8028 bytes).	0xAB 0xAA	No error.

00000040 ba 46 34 2b 21 00 00 20 04 00 00 00 c7 f4 20 8b  
00000050 71 c0 1c 9f 26 84 a7 9a b9 45 ce 8a 5f 33 f6 f1  
00002040 89 18 fd db c5 3b ef 49 11 35 bb 06 b5 25 23 53  
0000205b 7b e5 01 8a 74 4b 02 f6 98 d4 75 bf 03 2a c6 b0

Figure 5. Send page bytes 0x4C to 0x205B from the .msbl file.

0xAA 0x80 0x04 0x03 0x2A 0xC6 ... 0xA8 0xEE 0xFC* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0x205C to 0x406B from the .msbl file (8028 bytes).	0xAB 0xAA	No error.
--	--	-----------	-----------



0xAA 0x80 0x04 0x39 0x43 0x5B ... 0x12 0x3E 0x1B* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0x406C to 0x607B from the .msbl file (8028 bytes).	0xAB 0xAA	No error.
0xAA 0x80 0x04 0x39 0xAC 0xB8 ... 0x14 0x6F 0x13* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0x607C to 0x808B from the .msbl file (8028 bytes).	0xAB 0xAA	No error.
0xAA 0x80 0x04 0xAA 0xF3 0x95 ... 0x53 0x08 0xEA* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0x808C to 0xA09B from the .msbl file (8028 bytes).	0xAB 0xAA	No error.
0xAA 0x80 0x04 0xB1 0xBF 0x28 ... 0x89 0x86 0x35* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0xA09C to 0xC0AB from the .msbl file (8028 bytes).	0xAB 0xAA	No error.
0xAA 0x80 0x04 0x56 0x26 0x44 ... 0x72 0x1D 0x74* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0xC0AC to 0xE0BB from the .msbl file (8028 bytes).	0xAB 0xAA	No error.
0xAA 0x80 0x04 0x95 0x25 0x5E ... 0x43 0xF0 0x7F* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0xE0BC to 0x100CB from the .msbl file.	0xAB 0xAA	No error.
0xAA 0x80 0x04 0x94 0xEC 0xEC ... 0x0B 0x37 0x95* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0x100CC to 0x120DB from the .msbl file (8028 bytes).	0xAB 0xAA	No error.
0xAA 0x80 0x04 0xA6 0xEA 0x53 ... 0xC2 0x26 0x87* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0x120DC to 0x140EB from the .msbl file (8028 bytes).	0xAB 0xAA	No error.
0xAA 0x80 0x04 0x01 0x87	Bootloader flash. Send page bytes 0x140EC to 0x160FB	0xAB 0xAA	No error.

0x1D ... 0xBE 0xFF 0xEC* (CMD_DELAY = 680ms)	from the .msbl file (8028 bytes).		
0xAA 0x80 0x04 0xEA 0x65 0xD6 ... 0x63 0xCF 0x5F* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0x160FC to 0x1810B from the .msbl file (8028 bytes).	0xAB 0xAA	No error.
0xAA 0x80 0x04 0x4F 0x7E 0xEF ... 0x26 0x3D 0x47* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0x1810C to 0x1A11B from the .msbl file (8028 bytes).	0xAB 0xAA	No error.
0xAA 0x80 0x04 0xF0 0x64 0x3D ... 0xAA 0x65 0x47* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0x1A11C to 0x1C12B from the .msbl file (8028 bytes).	0xAB 0xAA	No error.
0xAA 0x80 0x04 0x54 0xDF 0xAB ... 0xA4 0xA9 0xE5* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0x1C12C to 0x1E13B from the .msbl file (8028 bytes).	0xAB 0xAA	No error.
0xAA 0x80 0x04 0xB0 0x0A 0xA5 ... 0xA6 0x24 0xA8* (CMD_DELAY = 680ms) Bootloader flash. Send page bytes 0x1E13C to 0x2014B from the .msbl file (8028 bytes).			
	0xAB 0xAA	No error.	
0xAA 0x80 0x04 0xE6 0x8E 0x35 ... 0x25 0x11 0xD4* (CMD_DELAY = 680ms)*	Bootloader flash. Send page bytes 0x2014C to 0x2215B from the .msbl file (8028 bytes).	0xAB 0xAA	No error.
0xAA 0x80 0x04 0xA0 0x88 0x34 ... 0x40 0x82 0x2E* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0x2215C to 0x2416B from the .msbl file (8028 bytes).	0xAB 0xAA	No error.
0xAA 0x80 0x04 0xD5 0x2F	Bootloader flash. Send page bytes 0x2416C to 0x2617B	0xAB 0xAA	No error.

0xAC ... 0xAC 0xCE 0x69* (CMD_DELAY = 680ms)	from the .msbl file (8028 bytes).		
0xAA 0x80 0x04 0x06 0x6C 0xF9 ... 0x9B 0xF2 0x91* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0x2617C to 0x2818B from the .msbl file (8028 bytes).	0xAB 0xAA	No error.
0xAA 0x80 0x04 01Bx 0xAD 0x09 ... 0xD3 0x7B 0xE6* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0x2818C to 0x2A19B from the .msbl file (8028 bytes).	0xAB 0xAA	No error.
0xAA 0x80 0x04 0xF3 0x08 0x6A ... 0x9C 0x46 0x31* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0x2A19C to 0x2C1AB from the .msbl file (8028 bytes).	0xAB 0xAA	No error.
0xAA 0x80 0x04 0x1E 0xCE 0x62 ... 0x99 0x60 0xC0* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0x2C1AC to 0x2E1BB from the .msbl file (8028 bytes).	0xAB 0xAA	No error.
0xAA 0x80 0x04 0xEF 0x61 0x8D ... 0xD9 0x40 0x92* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0x2E1BC to 0x301CB from the .msbl file (8028 bytes).	0xAB 0xAA	No error.
0xAA 0x80 0x04 0xA0 0x25 0x37 ... 0xA7 0x15 0xAA* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0x301CC to 0x321DB from the .msbl file (8028 bytes).	0xAB 0xAA	No error.
0xAA 0x80 0x04 0xF4 0xB0 0x84 ... 0x41 0x8B 0x6A* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0x321DC to 0x341EB from the .msbl file (8028 bytes).	0xAB 0xAA	No error.
0xAA 0x80 0x04 0x7D 0x02 0x6E ... 0x87 0x83 0x70*	Bootloader flash. Send page bytes 0x341EC to 0x361FB from the .msbl file (8028 bytes).		No error.

(CMD_DELAY = 680ms)			
0xAA 0x80 0x04 0x1E 0x8D 0x35 ... 0x18 0x19 0xD8* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0x361FC to 0x3820B from the .msbl file (8028 bytes).		No error.
0xAA 0x80 0x04 0xEE 0x60 0xBA ... 0x5C 0x27 0x13* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0x3820C to 0x3A21B from the .msbl file (8028 bytes).		No error.
0xAA 0x80 0x04 0x58 0xF0 0xE9 ... 0xAF 0x7E 0x79* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0x3A21C to 0x3C22B from the .msbl file (8028 bytes).		No error.
0xAA 0x80 0x04 0xD3 0x6F 0x26 ... 0x30 0xCB 0x5C* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0x3C22C to 0x3E23B from the .msbl file (8028 bytes).		No error.
0xAA 0x80 0x04 0x8F 0xAA 0xB6 ... 0xF6 0x66 0xEF* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0x3E23C to 0x4024B from the .msbl file (8028 bytes).		No error.
0xAA 0x80 0x04 0x8A 0x0B 0xD3 ... 0xE2 0x24 0x74* (CMD_DELAY = 680ms)	Bootloader flash. Send page bytes 0x4024C to 0x4225B from the .msbl file (8028 bytes).		No error.




Figure 6. Sequence to enter application mode.

Alternately, the MAX32674C can be commanded to application mode.<sup>†</sup>


0xAA 0x01 0x00 0x00 <sup>†</sup> (CMD_DELAY = 1.5s)	Set mode to 0x00 for application mode.	0xAB 0xAA	No error.
0xAA 0x02 0x00+	Read mode.	0xAB 0xAA 0x00	No errors. Mode is application.

\*Mandatory

†Recommended

For host micros that have limited RAM or limited I2C message sizes, the AlgoHub allows for the pages to be loaded with among several messages. The table below illustrates how to flash an .msbl file which has 5 pages using the partial page writing. In this example, the partial page size is set to 4000 bytes.

**Table 21. Annotated I2C Trace for Flashing the Application using the Partial Page Writing Command**

HOST COMMAND	COMMAND DESCRIPTION	READ MAX32674C RESPONSE	RESPONSE DESCRIPTION
<p>Sequence the MAX32674C to enter bootloader mode. *</p>  <p>Figure 7. Sequence to enter bootloader mode.</p>			
0xAA 0x02 0x00†	Read mode.	0xAB 0xAA 0x08	No error. Mode is bootloader.
0xAA 0x81 0x01†	Read bootloader page size.	0xAB 0xAA 0x20 0x00	No error. Page size is 8192.
0xAA 0x80 0x02 0x00 0x05*	Bootloader flash. Set the “number of pages” to 5 based on the value at byte 0x44 from the application .msbl file.	0xAB 0xAA	No error.
0xAA 0x80 0x00 [11 bytes]*	Bootloader flash. Set the initialization vector bytes to the 0x28 to 0x32 values from the .msbl file.	0xAB 0xAA	No error.
0xAA 0x80 0x01 [16 bytes]*	Bootloader flash. Set the authentication bytes to the 0x34 to 0x43 values from the .msbl file.	0xAB 0xAA	No error.
0x80 0x06 0x0F 0xA0*	Set partial page load size as 4000 (0x0FA0)	0xAB 0xAA	No error.
0xAA 0x80 0x03* (CMD_DELAY = 1400ms)	Bootloader flash. Erase application.	0xAB 0xAA	No error.
0xAA 0x80 0x04 [4000 bytes]*	Bootloader flash. Send first partial bytes of the first	0xAB 0xAA	No error.

(CMD_DELAY = 680ms)	page, 0x4C to 0xFEB from the .msbl file.		
0xAA 0x80 0x04 [4000 bytes]* (CMD_DELAY = 680ms)	Bootloader flash. Send second partial bytes of the first page, 0xFEC to 0x1F8B from the .msbl file.	0xAB 0xAA	No error.
0xAA 0x80 0x04 [208 bytes]* (CMD_DELAY = 680ms)	Bootloader flash. Send the remaining partial page bytes of the first page, 0x1F8C to 0x205B from the .msbl file	0xAB 0xAA	No error.
0xAA 0x80 0x04 [4000 bytes]* (CMD_DELAY = 680ms)	Bootloader flash. Send first partial bytes of the second page, 0x205C to 0x2FFB from the .msbl file.	0xAB 0xAA	No error.
0xAA 0x80 0x04 [4000 bytes]* (CMD_DELAY = 680ms)	Bootloader flash. Send second partial bytes of the second page, 0x2FFC to 0x3F9B from the .msbl file.	0xAB 0xAA	No error.
0xAA 0x80 0x04 [208 bytes]* (CMD_DELAY = 680ms)	Bootloader flash. Send the remaining partial page bytes of the second page, 0x3F9C to 0x406B from the .msbl file	0xAB 0xAA	No error.
0xAA 0x80 0x04 [4000 bytes]* (CMD_DELAY = 680ms)	Bootloader flash. Send first partial bytes of the third page, 0x406C to 0x500B from the .msbl file.	0xAB 0xAA	No error.
0xAA 0x80 0x04 [4000 bytes]* (CMD_DELAY = 680ms)	Bootloader flash. Send second partial bytes of the third page, 0x500C to 0x5FAB from the .msbl file.	0xAB 0xAA	No error.
0xAA 0x80 0x04 [208 bytes]* (CMD_DELAY = 680ms)	Bootloader flash. Send the remaining partial page bytes of the third page, 0x5FAC to 0x607B from the .msbl file	0xAB 0xAA	No error.

0xAA 0x80 0x04 [4000 bytes]* (CMD_DELAY = 680ms)	Bootloader flash. Send first partial bytes of the fourth page, 0x607C to 0x701B from the .msbl file.	0xAB 0xAA	No error.
0xAA 0x80 0x04 [4000 bytes]* (CMD_DELAY = 680ms)	Bootloader flash. Send second partial bytes of the fourth page, 0x701C to 0x7FBB from the .msbl file.	0xAB 0xAA	No error.
0xAA 0x80 0x04 [208 bytes]* (CMD_DELAY = 680ms)	Bootloader flash. Send the remaining partial page bytes of the fourth page, 0x7FBC to 0x808B from the .msbl file	0xAB 0xAA	No error.
0xAA 0x80 0x04 [4000 bytes]* (CMD_DELAY = 680ms)	Bootloader flash. Send first partial bytes of the last page, 0x808C to 0x902B from the .msbl file.	0xAB 0xAA	No error.
0xAA 0x80 0x04 [4000 bytes]* (CMD_DELAY = 680ms)	Bootloader flash. Send second partial bytes of the last page, 0x902C to 0x9FCB from the .msbl file.	0xAB 0xAA	No error.
0xAA 0x80 0x04 [208 bytes]* (CMD_DELAY = 680ms)	Bootloader flash. Send the remaining partial page bytes of the last page, 0x9FCC to 0xA09B from the .msbl file	0xAB 0xAA	No error.




Figure 8. Sequence to enter application mode.

Alternately, the MAX32674C can be commanded to application mode.<sup>†</sup>

0xAA 0x01 0x00 0x00 <sup>†</sup> (CMD_DELAY = 1.5s)	Set mode to 0x00 for application mode.	0xAB 0xAA	No error.
0xAA 0x02 0x00+	Read mode.	0xAB 0xAA 0x00	No errors. Mode is application.

\*Mandatory

<sup>†</sup>Recommended



## SpO<sub>2</sub> Coefficients for Final Factor

Due to variations in the physical design and optical cover lens of the final product, a calibration data collection procedure for SpO<sub>2</sub> is required to be performed once in a controlled environment. This procedure is important to ensure the quality of the SpO<sub>2</sub> calculation. This step is typically performed in a standard lab using the final form factor (with cover lens) with a reference SpO<sub>2</sub> device to determine three SpO<sub>2</sub> calibration coefficients: a, b, and c. The details of the SpO<sub>2</sub> calibration data collection and SpO<sub>2</sub> coefficient derivation procedure are described in the [Guidelines for SpO<sub>2</sub> Measurement Using the Maxim MAX32664 Sensor Hub](#) application note.

Once the three SpO<sub>2</sub> calibrations coefficients are obtained, they need to be loaded to the sensor hub every time prior to starting the algorithm using the “write SpO<sub>2</sub> calibration coefficients command.”

The SpO<sub>2</sub> calibrations coefficients need to be converted to a 32-bit integer format using the following:

- $A_{int32} = \text{round}(10^5 \times a)$
- $B_{int32} = \text{round}(10^5 \times b)$
- $C_{int32} = \text{round}(10^5 \times c)$

The SpO<sub>2</sub> calibration coefficients may be stored in the host flash separately and loaded to the sensor hub after every reset.

## APIs for Sleep, Shutdown

Summarized below are the commands and methods to place the MAX32674C, AFE, accelerometer into sleep or shutdown.

**Table 22. Sleep, Shutdown I2C Message Protocol Definitions**

COMMAND NAME	HOST COMMAND TO MAX32674C	DESCRIPTION
--------------	------------------------------	-------------

MAX32674C sleeps when idle and MFIO is low		MAX32674C uses deep-sleep for low-powered mode when idle and MFIO low.
MAX32674C shutdown	0xAA 0x01 0x00 0x01	Place the MAX32674C into shutdown. Restart by power cycling or pulsing RSTN.
MAX32674C hard reset	Use MFIO and RSTN pins according to section 2.	
WDT in MAX32674C		Not implemented.
MAX86176 AFE reset. Use AFE write register command. SensorHub.	0xAA 0x40 0x06 0x10 0x01	Write 0x01 to MAX86176 register 0x10 (System Configuration) to put the reset the MAX86176.
MAX86176 AFE shutdown. Use AFE write register command. SensorHub.	0xAA 0x40 0x06 0x10 0x02	Write 0x02 to MAX86176 register 0x10 (System Configuration) to put the MAX86176 into shutdown (SHDN) mode.
LIS2DS12 wake on activity or wake on double/single tap. Use accel write register command. SensorHub.	0xAA 0x40 0x04 [reg_addr] [value]	See AN4748 LIS2DS12 for register settings. Accel interrupt pin can be connected to the host GPIO.

## MAX32674C Processing Capabilities

The MAX32674C IC hardware is the same as the MAX32670.

- MIPS: Arm Cortex-M4 with FPU: 1.27 Dhrystone MIPS/MHz
- RAM: 160KB SRAM
- Flash: 384KB Flash Memory with Error Correction
- CPU Frequency: 100MHz

## Heart Rate Algorithm Performance

**Table 23. Heart Rate Algorithm Performance**

Wearable Heart Rate Monitoring – WHRM		
Category	Features	Specifications
Algorithm	Measurement principle:	Optical PPG signal from wrist, 3D Axis Accelerometer
	Measurement range:	HR: [30 240] BPM Cadence (steps per minute): [90 360]
	Measurement accuracy:	HR: Accuracy Definition -> within +/-10% error band vs. reference (chest strap) Resting: 94 – 100% Walking: 93 - 99% Biking: 91 - 97% Running: 91 - 97% Daily Life: 90 - 100% Step counting: Accuracy Definition: 100% – (Absolute Percent Error) Treadmill walking: 89 – 95% Treadmill running: 86 – 92% Outdoor walking: 80 – 90%

		Activity Classification: Accuracy Definition: 100% – (Absolute Percent Error) Rest: 87 - 93% Treadmill walking: 93 - 99% Treadmill running: 90 - 95% Outdoor walking: 91 - 97% Outdoor biking: 80 - 90% Energy consumption Kcal: calculated according to ACSM & ADA
	Reference measurement device:	HR: ECG based chest strap IBI: ECG based chest strap (sampling rate min 1kHz)
	Average response time	25Hz, first response time 15sec
	Inputs:	Single/Multiple Channel PPG signal 3-axis accelerometer signals
	Built-in features:	Activity Classifier Built-in Step Counter Motion compensation of PPG for accurate HR estimation Inter-beat interval estimator Energy expenditure estimation
Measurement Positions	Wrist, Ear, Finger, Chest, Abdomen	Sports and daily life activities
Sensor & Signal Requirements	LED requirements:	Please refer to “Reference Design Document” for details
	Perfusion index range:	Minimum AC to calculate HR is 20nA with average 0.8% PI
	Sampling rate:	25 Hz
Calibration	Calibration:	Algorithm activity classifier is tuned for the sensor placement on wrist, a calibration might be required to train algorithm to improve its performance for another body location.

## SpO2 Algorithm Performance

**Table 24. SpO2 Algorithm Performance**

SpO <sub>2</sub> on Wrist		
Category	Features	Specifications
Algorithm	Measurement principle:	Optical PPG signal from wrist
	Measurement range:	70 – 100% SpO <sub>2</sub>
	Measurement accuracy:	RMSE ≤ 3.5% as required by FDA for reflective mode pulse oximeters
	Clinical test:	Certified calibration lab
	Measurement time:	30 – 60 sec, on-demand operations, one shot
	Inputs:	Red and Infrared PPG signals 3-axis accelerometer signals
	Built-in features:	Precise motion detector Automatic AFE setting adjustment for optimum PPG quality Signal conditioning Signal selection according to signal quality for discarding noisy signal portions Adjustable time-out duration Adjustable confidence threshold
Measurement Positions	Standing:	Arm is kept horizontal at the level of heart; palm is facing the floor
	Sitting	Arms are placed on a table Arms are crossed, arm with wristwatch is above the other arm
	Lying down:	Arms are horizontal

Sensor & Signal Requirements	LED requirements:	Center wavelength shift $\leq \pm 5\text{nm}$ LED full width at half maximum (FWHM) $\leq 20\text{nm}$
	Perfusion index range:	PI $\geq 0.05\%$
	Sampling rate:	25 Hz
	Optical layout design:	Please refer to "Design Guide for SpO2 Measurement"
Calibration	Calibration lab:	Lab calibration is required for wearable's finished industrial design. Please refer to lab calibration procedure guide "Design Guide for SpO2 Measurement"

## References

[MAXREFDES104#: Health Sensor Platform 3.0 \(Design Files, Firmware, Software, User Guide, Sample Host Code\)](#)

[Validation and Performance of a Wearable Heart rate Monitoring Algorithm](#)

[Guidelines For Spo2 Measurement](#)

[Search Knowledge Base for FAQ](#)

[Technical Support](#)